

### 3. SPECIFICITY

The antiserum used in the assay is highly specific for the detection of 17 $\alpha$ -hydroxyprogesterone. The following naturally occurring steroids were screened for cross-reactivity.

Compound	Concentration (ng/ml)	% Cross Reactivity
Estrone	10	0.10
Estradiol	40	0.02
Estriol	400	0.02
Progesterone	400	0.30
Cortisol	8000	0.01
Testosterone	200	0.01
Androstenedione	1000	0.01
DHEA-S	100000	0.00

2020-03-09



## NEONATAL 17-OH PROGESTERONE ELISA

Catalog No.: PG709S (96 Tests)

### INTENDED USE

The Neonatal 17- hydroxyprogesterone (17-OHP) ELISA kit is intended for the quantitative determination of 17-OHP in Human whole blood.

### PRINCIPLE OF THE TEST

The Neonatal 17 $\alpha$ -OH Progesterone (N-17-OHP) is a competitive solid phase ELISA in blood spot dried on WHATMAN type 903 filter paper. In the assay, goat anti-rabbit IgG-coated wells are incubated with blood spots, 50 $\mu$ l N-17-OHP-HRP conjugate, and 50 $\mu$ l rabbit anti-N-17-OHP. HRP-labeled-N-17-OHP competes with N-17-OHP in the standards and sample for a fixed number of binding sites of the specific anti-17-OHP antibody. Thus, the amount of 17-OHP peroxidase conjugate immunologically bound to the well progressively decreases as the concentration of 17-OHP in the specimen increases. Unbound 17-OHP is then removed by the washing steps. TMB Substrate is added, resulting in the development of blue color. The color development is stopped with the addition of stop solution, and the absorbance is spectrophotometrically measured at 450nm. A standard curve is prepared relating color intensity to the concentration of 17-OHP.

MATERIALS PROVIDED	96 Tests
1. Microwells coated with Goat anti-Rabbit IgG	12x8x1
2. 17-OH Progesterone Standards, 6 Dried Blood Spots	6 DBS
3. 17-OH Progesterone Controls, 2 Dried Blood Spot per level, per card	6 DBS
4. Enzyme Conjugate, 20X: 1 vial	0.35 mL
5. 17-OH Progesterone Assay Diluent: 1 bottle	7 mL
6. Rabbit anti-17-OH Progesterone Reagent: 1 bottle (ready to use)	7 mL
7. TMB Substrate: 1 bottle (ready to use)	12 mL
8. Stop Solution: 1 bottle (ready to use)	12 mL
9. Wash Concentrate 20X: 1 bottle	25 mL


### MATERIALS NOT PROVIDED

- 1/8" inch hole punch
- Nutator capable of fixed speed rotation
- Distilled or deionized water
- Precision pipettes
- Disposable pipette tips
- ELISA reader capable of reading absorbance at 450 nm
- Absorbance paper or paper towel
- Graph paper

### STORAGE AND STABILITY

- Store the kit at 2-8 $\circ$ C.
- Keep microwells sealed in a dry bag with desiccants.
- The reagents are stable until expiration of the kit.
- Do not expose test reagents to heat, sun or strong light.

Cat#: PG709S (96 Tests)



CEpartner4U

3951 DB; 13 NL  
Tel +31 (0) 6-516.536.26

**WARNINGS AND PRECAUTIONS**

- Potential biohazardous materials: the calibrator and controls contain human source components, which have been tested and found non-reactive for hepatitis B surface antigen as well as HIV antibody with FDA licensed reagents. However, there is no test method that can offer complete assurance that HIV, Hepatitis B virus or other infectious agents are absent. These reagents should be handled at the Biosafety Level 2, as recommended in the Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories." 1984.
- Optimal results will be obtained by strict adherence to the test protocol. Precise pipetting as well as following the exact time and temperature requirements is essential.
- Do not pipette by mouth. Do not smoke, eat, or drink in the areas in which specimens or kit reagents are handled.
- The components in this kit are intended for use as an integral unit. The components of different lots should not be mixed.

**SPECIMEN COLLECTION AND HANDLING**

Follow the guidelines in the NCCLS publication LA4T7 for collecting blood samples in the neonatal screening program, copies of which can be obtained from: NCCLS, 771 E. Lancaster Ave, Villanova, PA 19085. Use WHATMAN type 903 filter paper. For samples screening for CAH, collect samples 3 to 5 days after birth. Use disposable lancets with tips less than 2.5 mm to prick the medial or lateral sides of the bottom of the heel. Allow a drop of blood to form with sufficient volume to fill a 5/8-inch diameter spot on filter paper. Gently touch the drop of blood with the filter paper. **DO NOT PRESS AGAINST THE SKIN. DO NOT TOUCH SPOTTED AREA.** Suspend spotted papers horizontally and allow drying at room temperature for a minimum of 3 hours. Avoid spots touching other surfaces and keep away from direct light. The samples should be transported to the laboratory within 24 hours after collection in appropriate storage container. The laboratory should store the specimens at 2-8 °C protected from moisture and direct light.

The dried blood spots are stable for at least 3 weeks at 2-8 °C protected from light and moisture. Reject samples with the following conditions:

- Specimens not collected on WHATMAN type 903 filter paper.
- Blood spots not completely saturated on both sides.
- Blood spots with appearance of caking or clotting.
- Blood spots with appearance of moisture

**REAGENT PREPARATION**

- Working Enzyme Conjugate: Prepare 1X working dilution at 1:20 with assay diluents as needed, e.g. 0.1ml of the stock conjugate in 1.9ml of assay diluent for 40 wells. The diluted conjugate has to be used on the same day.
- Wash Buffer: Prepare 1X Wash Buffer by adding the contents of the bottle (25ml, 20X) to 475ml of distilled water. Store at room temperature (20-25°C).

**ASSAY PROCEDURE**

- Place the desired number of coated strips into the holder.
- Punch out 1/8" blood spot out of each standard and samples into the assigned wells. (NOTE: Do not punch blood spots from areas that are printed or that are near the edge of the blood spot).
- Add 50µl of working dilution N-17-OHP enzyme conjugate to all the wells.
- Add 50 µl of N-17-OHP Antibody Reagent to all wells.
- Shake the microplate gently for 20-30 seconds to mix. Make sure that all blood spots are fully submerged in the liquid and not stuck to the walls of the wells.
- Cover the microplate and put it on the nutator for overnight incubation (18-20 hours). Alternatively, microplate can be incubated for 4hours, on a plate shaker @1500rpm, with an intermittent 5 seconds agitation, for every 30 seconds break.
- Remove the contents of the wells by decantation or aspiration. Make sure all the blood dots are removed at this point. Rinse the wells 5 times with 1X wash buffer. Strike the wells sharply to remove residual wash buffer droplets.
- Add 100µl of TMB substrate to each well.

- Cover the microplate and incubate for 15 minutes at room temperature.
- Add 50µl of stop solution to each well and gently mix until a uniform color, in each well, is obtained.
- Read the absorbance in each well at 450nm within 15 minutes after adding the stop solution.

**CALCULATION OF RESULTS**

A standard curve is constructed as follows:

- Calculate the average absorbance values for each set of standards and patient samples
- To construct the standard curve, plot the mean absorbance of each N-17-OHP standards (vertical axis) against its concentration in ng/ml (horizontal axis)
- Draw the best-fit curve through the plotted points.
- Read the absorbance for each unknown sample from the curve to determine the corresponding concentration of N-17-OHP.

**Example of a Typical Standard Curve**

	OD450nm	Conc. (ng/mL)
Std 1	2.399	0
Std 2	2.080	5
Std 3	1.372	14
Std 4	0.698	35
Std 5	0.405	70
Std 6	0.235	140

**EXPECTED VALUES**

We recommend each laboratory to establish its own normal ranges, for the population it serves. Until then, literature values may be used as guidelines.

**REPORTABLE RANGE:** Analytical Range = 5 - 140 ng/ml Samples that fall within the calibration curve should be reported as such. Samples that fall outside the calibration curve should be reported <5 ng/mL or >140) ng/ml.

**EXPECTED VALUES**

We recommend each laboratory to establish its own normal ranges, for the population it serves. Until then, literature values may be used as guidelines.

- Sensitivity**

The sensitivity was determined by calculating the mean plus 2SD of the standard zero point tested 24 times in the same run.

Serum	No. of Replicates	Mean (ng/mL)	Standard Deviation	Mean + 2SD (Sensitivity)
Zero Standard	24	0.229	0.299	0.827

- Precision**

**Intra-Assay**

Serum	No. of Replicates	Mean (ng/mL)	Standard Deviation	Coefficient of Variation (%)
1	16	26.8	2.398	8.96%
2	16	43.3	2.408	5.56%
3	16	90.6	4.145	4.58%

**Inter-assay**

Serum	No. of Replicates	Mean (ng/ml)	Standard Deviation	Coefficient of Variation (%)
1	20	28.5	2.512	8.83%
2	20	44.0	2.736	6.21%
3	20	89.0	4.479	5.03%