DCM051-10
Ed. 03/2024

CEA ELISA

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CEA ELISA

for routine analysis

Direct immunoenzymatic determination of Carcinoembryonic Antigen (CEA) in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO051

INTENDED USE

Dia.Metra CEA ELISA kit is a direct solid phase enzyme immunoassay for quantitative determination of Carcinoembryonic Antigen (CEA) in human serum or plasma.

CEA ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Carcinoembryonic antigen (CEA) is a glycoprotein, with a molecular weight of 180 kD, involved in cell adhesion. It is normally produced during fetal development, but the production of CEA stops before birth. Therefore, it is not usually present in the blood of healthy adults, although levels are raised in heavy smokers. CEA was identified in human colon cancer tissue extracts. It was later found that serum from individuals with colorectal and other carcinomas had higher levels of CEA than healthy individuals and can be used to monitor the response to colon cancer treatment.

CEA and related genes make up the CEA family belonging to the immunoglobulin superfamily.


The most frequent cancer which causes an increased CEA is cancer of the colon and rectum. Others include cancers of the pancreas, stomach, breast, lung, and certain types of thyroid and ovarian cancer. Benign conditions which can elevate CEA include smoking, infections, inflammatory bowel disease, pancreatitis, cirrhosis of the liver, and some benign tumors in the same organs in which an elevated CEA indicates cancer. Chemotherapy and radiation therapy can cause a temporary rise in CEA due to the death of tumor cells and release of CEA into the blood stream. Benign tumours does not usually cause an increase above 10 ng/mL.

2. PRINCIPLE

In this method, the calibrators, the patient specimens and/or the controls (containing the native CEA antigen) are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labelled antibodies are then added and the reactants mixed: these antibodies have high affinity and specificity and are directed against distinct and different epitopes of CEA.

Reaction between the various CEA antibodies and native CEA occurs in the microwells without competition or steric hindrance, forming a soluble sandwich complex.

The interaction is illustrated by the following equation:

BtAb(m) = Biotinylated Monoclonal Antibody (Excess Quantity)

AgCEA = Native Antigen (Variable Quantity)

EAb = Enzyme labelled Antibody (Excess Quantity)

EAb-AgCEA-BtAb(m) = Antigen-Antibodies Sandwich Complex

ka = Rate Constant of Association

k-a = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

Streptavidin C.W. = Streptavidin immobilized on well
Immobilized complex = sandwich complex bound to the well

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by a washing step.

The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantified by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

ProClin <0.0015%

CAL0	REF	DCE002/5106-0
CAL1	REF	DCE002/5107-0
CAL2	REF	DCE002/5108-0
CAL3	REF	DCE002/5109-0
CAL4	REF	DCE002/5110-0
CAL5	REF	DCE002/5111-0

2. Conjugate (1 vial, 13 mL)

Antibodies anti CEA conjugated with Horseradish peroxidase (HRP) and anti CEA biotinylated, ProClin <0.0015% REF DCE002/5102-0

3. Coated Microplate (1 microplate breakable)

Microplate coated with streptavidin REF DCE002/5103-0

4. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB 0.26 g/L (avoid any skin contact), ProClin <0.0015% REF DCE004-0

5. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact) REF DCE005-0

6. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L; Tween-20 55 g/L, ProClin >0.0015% REF DCE006-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser

Microplates reader (450 nm, 620-630 nm)

Note

Store all reagents at 2÷8°C in the dark.

Open the bag of reagent 3 (Coated Microplate) only when it is at room temperature and close immediately after use; once opened, it is stable until the expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators should be handled in the same manner as potentially infectious material.

- Some reagents (wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- **Classification according to Regulation (EC) No. 1272/2008 [CLP]**

Skin sensitivity, Category 1

Contains: ProClin 300

Warning

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of CEA from 5.0 to 250 ng/mL.
- CEA has a low clinical sensitivity and specificity as a tumour marker. Clinically an elevated CEA value alone is not of diagnostic value as a test for cancer and should only be used in conjunction with other clinical manifestations (observations) and diagnostic parameters. There are patients with colorectal cancer that do not exhibit elevated CEA values and elevated CEA values do not always change with progression or regression of disease. Smokers demonstrate a higher range of baseline values than non-smokers.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems, it is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.

Reagent	Calibrator	Sample	Blank
Calibrator C ₀ -C ₅	25 µL		
Samples		25 µL	
Conjugate	100 µL	100 µL	
Incubate at room temperature (22÷28°C) for 1 hours. Remove the contents from each well. Wash the wells 3 times with 300 µL of diluted Wash Solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₅)

The Calibrators are ready to use, are calibrated using a reference preparation which was assayed against the 1st IRP 73/601 and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	5	10	25	50	250

Once opened, the Calibrators are stable 6 months at 2-8°C.

6.2. Preparation of Wash Solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8°C.

6.3. Preparation of the Sample

The specimens shall be human serum or plasma; the usual precautions in the collection of venipuncture samples should be observed.

For accurate comparison to established normal values, a fasting morning serum sample should be obtained.

For serum preparation, the blood should be collected in a venipuncture tube without additives or anti-coagulants; then allow the blood to clot and centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2÷8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they may be stored at temperatures of -20°C for up to 30 days. Avoid repetitive freezing and thawing.

When assayed in duplicate, 0.050 mL of the specimen is required.

For sample with concentration over 250 ng/mL dilute the sample with Calibrator 0.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of CEA for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Note

The optical densities (O.D.s) of some calibrators and samples may be higher than 2.0, in such a case, they could be out of the measurement range of the microplate reader. It is therefore necessary, for O.D.s higher than 2.0, to perform a reading at 405 nm (=wavelength of peak shoulder) in addition to 450 nm (peak wavelength) and 620 (reference filter for the subtraction of interferences due to the plastic).

For microplate readers unable to read the plate at 3 wavelengths at the same time, it is advisable to proceed as follows:

- Read the microplate at 450 nm and at 620 nm.
- Read again the plate at 405 nm and 620 nm.
- Find out the wells whose ODs at 450 nm are higher than 2.0
- Select the corresponding ODs read at 405 nm and multiply these values at 405 nm by the conversion factor 3.0 (where $OD_{450}/OD_{405} = 3.0$), that is:
 $OD_{450\text{ nm}} = OD_{405\text{ nm}} \times 3.0$.

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

The OD of calibrator 5 should be ≥ 1.3 .

8.2. Mean Absorbance

Calculate the mean of the absorbance (E_m) for each point of the calibration curve (C_0 - C_5) and of each sample.

8.3. Calibration curve

Plot the mean value of absorbance (E_m) of the Calibrators (C_0 - C_5) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

8.4. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in ng/mL.

9. REFERENCE VALUES

Nearly 99% of non-smokers have CEA concentrations less than 5 ng/ml. Similarly 99% of smokers have concentrations less than 10 ng/ml (4).

	ng/mL
Non-smoker	<5 ng/mL
Smoker	<10 ng/mL

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (16x) of three different control sera in one assay. The within assay variability is $\leq 5.8\%$.

10.1.2. Inter Assay Variation

Between run variations was determined by replicate measurements (20x) of three different control sera in different lots of kit. The between assay variability is $\leq 9.9\%$.

10.2. Accuracy

The recovery of amounts of 10 - 20 - 40 ng/mL of CEA added to samples gave an average value ($\pm$ SD) of $100.7\% \pm 7.2\%$ with reference to the original concentrations.

10.3. Sensitivity

The lowest detectable concentration of CEA that can be distinguished from the Calibrator 0 is 0.11 ng/mL at the 95% confidence limit.

10.4. Specificity

In order to assess the specificity of the antibody pair used for the CEA Elisa assay, massive doses of related analytes were spiked in a pool of patient sera:

Analyte	Concentration	Cross Reaction
CEA		100 %
β -HCG	250 ng/mL	ND
AFP	5000 ng/mL	ND
CA-125	1000 U/mL	ND
PSA	1000 ng/mL	ND
CA 19-9	1000 U/mL	ND

10.5. Correlation

Dia.Metra CEA ELISA was compared to another commercially available CEA assay. 167 serum samples were analysed according in both test systems.

The linear regression curve was calculated:

$$y = 1.019x - 0.200$$

$$r^2 = 0.992$$

y = CEA Predicate kit

x = CEA Dia.Metra kit

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

- 9) Staab HJ, et al Am. J.Surgery 136:322-327 (1978)
- 10) Thomas P, et al Biochem Biophys Acta 1032:177-189 (1990)
- 11) Yamashita K, et al Cancer Research 47:3451-3459 (1987)
- 12) Hammerstrom S, et al Cancer Research 49:4852-58 (1989)
- 13) Ann.Inter.Med.1981;94:407-409

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