DCM016-15
Ed. 09/2024

CIC C1q ELISA

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CIC C1q ELISA

for routine analysis

Direct immunoenzymatic determination of CIC C1q in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO016

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use
For Laboratory Professional Use

CIC C1q ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of CIC C1q in human serum or plasma from an adult population.

2. CLINICAL SIGNIFICANCE

The complement system is a biochemical cascade of the immune system that helps clear pathogens from an organism. It is derived from many small plasma proteins that work together to form the primary end result of cytolysis by disrupting the target cell's plasma membrane.

Activation of this system leads to cytolysis, chemotaxis, opsonization, immune clearance, and inflammation, as well as the marking of pathogens for phagocytosis. The complement system consists of more than 35 soluble and cell-bound proteins, 12 of which are directly involved in the complement pathways. The proteins account for 5% of the serum globulin fraction. The complement proteins are synthesized mainly by hepatocytes; however, significant amounts are also produced by monocytes, macrophages, and epithelial cells in the gastrointestinal and genitourinary tracts.

C1q is involved in the classical complement pathway. The classical pathway is triggered by activation of the C1-complex (which consists of one molecule C1q and two molecules C1r and C1s), either by C1q's binding to antibodies from classes M and G, complexed with antigens, or by its binding C1q to the surface of the pathogen.

The complement system might play a role in many diseases with an immune component, such as Barraquer-Simons Syndrome Alzheimer's disease, asthma, lupus erythematosus, various forms of arthritis, autoimmune heart disease and multiple sclerosis. Deficiencies in the terminal pathway predispose to both autoimmune disease and infections (particularly meningitis).

There are many tests for the determination of CIC, included the test of precipitation with PEG, radial immunodiffusion, and cellular tests like the test of Ray cell. Does not exist one

procedure to determinate all types of immunocomplex; in commerce exist some test to determinate fragments of the complex (e.g. C1q and C3d) that have an important diagnostic function.

3. PRINCIPLE OF THE METHOD

CIC C1q ELISA kit is based on the binding of C1q-linked immunocomplexes to C1q adsorbed on microplate. In the first step, the samples are added to the microplate adsorbed with C1q; during the following incubation, C1q-fixing circulating immune complexes (CIC) bind to the C1q immobilized on the microplate. The microplate is washed for remove the unbound serum proteins.

In the second step, the anti-human IgG conjugated with horseradish peroxidase (HRP) is added; it binds to the immunocomplex fixed on the microplate. The washing step removes the unbound conjugate.

In the third step, the TMB Substrate is added, and this reacts with the conjugate fixed on the microplate, developing a colorimetric reaction. The quantity of CIC IgG complex is proportional to the colour intensity read at 450 nm wavelengths.

The immunocomplex concentration in the sample is calculated through a calibration cuve. Heat aggregate human gamma globulin per mL ($\mu\text{Eq/mL}$) is the unit of measure of the results.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (3 vials, 1.5 mL each)

ProClin >0.0015%

CAL0

REF DCE002/1606-0

CAL1

REF DCE002/1607-0

CAL2

REF DCE002/1608-0

2. Controls (2 vials, 1.5 mL each, ready to use)

74 mM Phosphate buffer, pH 7.4, BSA 0.5%, ProClin >0.0015%

Negative Control

REF DCE045/1601-0

Positive Control

REF DCE045/1602-0

3. Incubation Buffer (1 vial, 50 mL)

74 mM Phosphate buffer, pH 7.4. ProClin >0.0015%

REF DCE008-0

4. Conjugate (1 vial, 0.5 mL)
Anti human IgG conjugated with horseradish peroxidase (HRP). ProClin >0.0015%, BSA 0.1%

REF DCE002/1602-0

5. Conjugate Buffer (1 vial, 20 mL)
74 mM Phosphate buffer, pH 7.4. ProClin >0.0015%

REF DCE009-0

6. Coated Microplate (1 breakable microplate)
Microplate coated with C1q **REF DCE002/1603-0**

7. TMB Substrate (1 vial, 15 mL)
H₂O₂-TMB 0.26 g/L (*avoid any skin contact*).
ProClin <0.0015% **REF DCE004-0**

8. Stop Solution (1 vial, 15 mL)
Sulphuric acid 0.15 mol/L (*avoid any skin contact*)
REF DCE005-0

9. 10X Conc. Wash Solution (2 vials, 50 mL each)
0.2M Phosphate buffer, pH 7.4. ProClin >0.0015%
REF DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser
Precision Pipetting Devices
Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  All human source material used in the preparation of the reagents have been tested and found negative for antibodies to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, controls, incubation buffer, conjugate, conjugate buffer and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinisation and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analysers which are not equipped with disposable tips.**
For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months*.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

* in use stability data supports reagent stability when used three times within this period.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C (serum & plasma)	24 hours
Freeze/thaw cycles (serum & plasma)	1 cycle
-20°C (serum & plasma)	6 months

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer.

The calibrators are ready for use and have the following concentration:

	C ₀	C ₁	C ₂
µgEq/mL	0	16	64

The controls are ready to use; the concentration is printed on the label.

9.2. Preparation of the Conjugate

Dilute the Conjugate (reagent 4) 1:100 with Conjugate buffer (reagent 5). The exact quantity is proportional to the number of the assays.

Mix well and avoid formation of foam. Stable for 3 hours at room temperature (22 – 28°C).

9.3. Preparation of the Wash Solution

Dilute the content of the vial “10X Conc. Wash Solution” with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of C1q can be performed in human serum or samples.

All serum and plasma samples must be pre-diluted with incubation buffer (reagent 3):

Sample	10 µL
Incubation buffer (reagent 3)	500 µL

Mix gently; avoid using vortex.

9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₂), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Calibrator C ₀ -C ₂	100 µL		
Control		100 µL	
Diluted sample		100 µL	

Incubate 30 minutes at +37°C ±0.5°C.
Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

Diluted conjugate	100 µL	100 µL	
-------------------	--------	--------	--

Incubate 30 minutes at +37°C ±0.5°C.
Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Washing: follow the same indications of the previous point.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 15 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake gently the microplate.
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁶.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A best fit curve, **including Calibrator 0 is required**. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

12. MEASURING RANGE

The assay measuring range (AMR) is 3.9 – 36.2 µgEq/mL. Any value that reads below 3.9 µgEq/mL should be reported as “< 3.9 µgEq/mL”. Any value that reads above 36.2 µgEq/mL should be reported as “> 36.2 µgEq/mL”.

13. UNCERTAINTY OF MEASUREMENT

The uncertainty associated with the measurement of the assay calibrators (C₀ – C₂) was determined following CLSI EP29-A “Expression of Measurement Uncertainty in Laboratory Medicine; Approved Guideline”. The uncertainty for each calibrator level was assessed over a minimum of 5 days (1 run per day) by 3 different operators. Result from one representative lot:

Calibrator	Mean conc. (µgEq/mL)	Calibrator uncertainty U _c (µgEq/mL)
CAL0	< LoB	N/A
CAL1	16	2.0
CAL2	64	9.2

14. INTERPRETATION OF RESULTS

µgEq/mL of aggregates IgG	Interpretation
< 16	The sample should be considered negative
16 – 18	The sample should be graded equivocal and repeat testing / sampling should be performed according to internal practices
>18	The sample should be considered positive

The above index and interpretation should be considered as guidelines only; it is recommended that each laboratory establish its own index and interpretation based upon its own patient population.

15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1. Sensitivity and Specificity

The sensitivity and specificity were determined with guidance from CLSI EP24-A2 “Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves” using 103 negative and 8 positive samples run on a single reagent lot.

		DKO016		Total
		Positive	Negative	
Reference ELISA	Positive	7	1	8
	Negative	3	100	103
Total		10	101	111

Sensitivity: 88%

Specificity: 97%

15.2. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A2, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.4 µgEq/mL
Limit of Detection (LoD)	1.7 µgEq/mL
Limit of Quantitation (LoQ)	3.9 µgEq/mL

15.3. Trueness

Since no certified reference material is available, trueness has been demonstrated through assessment of sensitivity and specificity – refer to section 15.1.

15.4. Precision

Precision of the CIC C1q ELISA was determined by performing a complex precision study.

Repeatability: A total of 5 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from 1 representative lot is shown below:

Sample	n	Mean Conc. (µgEq/mL)	Repeatability	
			SD	CV%
1	75	6.1	0.9	14%
2	75	11.2	1.5	13%
3	75	22.2	2.1	10%
4	75	31.5	2.4	8%
5	75	54.4	5.6	10%

Reproducibility: A total of 5 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from 1 representative lot is shown below:

Sample	n	Mean Conc. (µgEq/mL)	Reproducibility	
			SD	CV%
1	75	6.1	1.1	18%
2	75	11.2	2.1	18%
3	75	22.2	3.3	15%
4	75	31.5	4.7	15%
5	75	54.4	6.7	12%

15.5. Linearity

Linearity was evaluated based on CLSI EP-06, “Evaluation of the Linearity of Quantitative Measurement Procedures”. A serial dilution of high and low serum samples was prepared.

For C1q concentration by CIC C1q ELISA, the measurement procedure shows linearity for the interval from 3.9 to 36.2 µgEq/mL within the allowable deviation of linearity (ADL) of ±15 %.

15.6. Interference

The following substances do not interfere with a bias of >15% in the CIC C1q ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	8.3 g/dL
Triglyceride	500 mg/dL

15.7. Serum-plasma study

The CIC C1q ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP35-Ed1 "Assessment of Equivalence or Suitability of Specimen Types for Medical Laboratory Measurement Procedures". A total of 20 samples to cover the assay range were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (µEq/mL) [95% CI]	Correlation coefficient (r)
SST	1.09 [0.99 - 1.18]	-0.5 [-1.74 - 0.78]	0.99
Lithium Heparin	1.07 [0.96 - 1.17]	-1.0 [-2.40 - 0.45]	0.98
Sodium Heparin	1.02 [0.88 - 1.17]	-0.4 [-2.37 - 1.64]	0.96
EDTA	1.09 [1.01 - 1.17]	-0.8 [-1.86 - 0.35]	0.99

16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁷. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

18. BIBLIOGRAPHY

1. Triolo G., et al J.Clin.Lab.Immunol 13, 35-39 (1984)
2. Rong-jia Xu et al J.Immunol.Met. 135, 225-231(1990)
3. Menzel J.E.,et al J.Immunol. Met. 138, 16 (1991)
4. Muso E, et al Nippon Jinzo Gakkai Shi.36(4):345-54 (1994)
5. Yoshinoya S,et al J Clin Lab Immunol. 38(4):161-73 (1992)
6. Basic QC Practices On-line Course; <http://www.Westgard.com>.
7. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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