## Instructions for Use

# CMV IgM ELISA

IVD

CE 0197

REF EIA-3469

$\Sigma$  96

**Please use only the valid version of the Instructions for Use provided with the kit.  
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.  
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

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1	INTRODUCTION .....	2
2	PRINCIPLE OF THE TEST .....	2
3	WARNINGS AND PRECAUTIONS .....	3
4	REAGENTS.....	4
5	SPECIMEN COLLECTION AND PREPARATION.....	5
6	ASSAY PROCEDURE .....	6
7	CALCULATION OF RESULTS .....	8
8	QUALITY CONTROL.....	9
9	PERFORMANCE CHARACTERISTICS.....	9
10	LIMITATIONS OF USE.....	11
11	LEGAL ASPECTS .....	11
12	REFERENCES / LITERATURE.....	11

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## 1 INTRODUCTION

### 1.1 Intended Use

The **DRG Cytomegalie Virus (CMV) IgM Enzyme Immunoassay Kit** provides materials for the **quantitative** and **qualitative** determination of IgM-class antibodies to Cytomegalie Virus (CMV) in human serum plasma (EDTA, heparin or citrate plasma).

**This assay is intended for in vitro diagnostic use only.**

### 1.2 Summary and Explanation

Cytomegalovirus (CMV) is a member of the herpesvirus group (Betasubfamily, DNA virus of 150-200 nm). These viruses share a characteristic ability to remain dormant within the body over a long period. Initial CMV infection, which may have few symptoms, is always followed by a prolonged, inapparent infection during which the virus resides in cells without causing detectable damage or clinical illness. Severe impairment of the body's immune system by medication or disease consistently reactivates the virus from the latent or dormant state.

CMV is found universally throughout all geographic locations and socioeconomic groups, and infects between 50% and 85% of adults.

CMV infection is more widespread in developing countries and in areas of lower socioeconomic conditions.

For the vast majority of people, CMV infection is not a serious problem, but it is to certain high-risk groups: the unborn baby during pregnancy, people who work with children, and

immunocompromised persons, such as organ transplant recipients and persons infected with HIV.

The presence of virus resp. infection may be identified by Microscopy, PCR, Serology: CBR and detection of antibodies by ELISA.

IgM antibodies are the first to be produced by the body in response to a CMV infection. They are present in most individuals within a week or two after the initial exposure. IgM antibody production rises for a short time period and then declines. After several months, the level of CMV IgM antibody usually falls below detectable levels. Additional IgM antibodies are produced when latent CMV is reactivated.

IgG antibodies are produced by the body several weeks after the initial CMV infection and provide protection from primary infections. Levels of IgG rise during the active infection, then stabilize as the CMV infection resolves and the virus becomes inactive. After a person has been exposed to CMV, he or she will have some measurable amount of CMV IgM antibody in their blood for the rest of their life. CMV IgM antibody testing can be used, along with IgM testing, to help confirm the presence of a recent or previous CMV infection.

## 2 PRINCIPLE OF THE TEST

The **DRG Cytomegalie Virus (CMV) IgM ELISA Kit** is a solid phase enzyme-linked immunosorbent assay (ELISA). Patient samples are diluted with *Sample Diluent* and additionally incubated with *IgG-RF-Sorbent* to eliminate competitive inhibition from specific IgG. This pretreatment avoids false negative results.

Microtiter wells as a solid phase are coated with inactivated grade 2 Cytomegalovirus (CMV) antigen (strain AD-169).

**Diluted patient** specimens and **ready-for-use controls** are pipetted into these wells. During incubation Cytomegalie Virus (CMV)-specific antibodies of positive specimens and controls are bound to the immobilized antigens.

After a washing step to remove unbound sample and control material horseradish peroxidase conjugated anti-human IgM antibodies are dispensed into the wells. During a second incubation this anti-IgM conjugate binds specifically to IgM antibodies resulting in the formation of enzyme-linked immune complexes.

After a second washing step to remove unbound conjugate the immune complexes formed (in case of positive results) are detected by incubation with TMB substrate and development of a blue color. The blue color turns into yellow by stopping the enzymatic indicator reaction with sulfuric acid.

The intensity of this color is directly proportional to the amount of Cytomegalie Virus (CMV)-specific IgM antibody in the patient specimen. Absorbance at 450 nm is read using an ELISA microtiter plate reader.

### 3 WARNINGS AND PRECAUTIONS

- This kit is for in vitro diagnostic use only. For professional use only.
- Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Avoid contact with Stop Solution containing 0.2 mol/L H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
- TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided
- Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
- Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
- Never pipette by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according the national biohazard safety guideline or regulation.
- For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

## 4 REAGENTS

### 4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;  
Wells coated with inactivated grade 2 Cytomegalovirus (CMV) antigen (strain AD-169).  
(incl. 1 strip holder and 1 cover foil)
  2. **Sample Diluent** \*, 1 vial, 100 mL, ready to use,  
colored yellow; pH  $7.2 \pm 0.2$ .
  3. **IgG-RF-Sorbent**\*, 1 vial, 6.5 mL, ready to use,  
colored yellow;  
Contains anti-human IgG-class antibody.
  4. **Standard (Standard 1-3)**\*, 3 vials, ready to use;  
Standard 1: 2.0 mL; Standard 2-3: 1 mL each  
Concentrations: **50, 200, 400 DU/mL**,  
colored yellow, white caps.
  5. **Pos. Control** \*, 1 vial, 1.0 mL, ready to use;  
colored yellow, red cap.
  6. **Neg. Control** \*, 1 vial, 2.0 mL, ready to use;  
colored yellow, yellow cap.
  7. **Enzyme Conjugate** \*, 1 vial, 20 mL, ready to use,  
colored red,  
antibody to human IgM conjugated to horseradish peroxidase.
  8. **Substrate Solution**, 1 vial, 14 mL, ready to use,  
Tetramethylbenzidine (TMB).
  9. **Stop Solution**, 1 vial, 14 mL, ready to use,  
contains 0.2 mol/L H<sub>2</sub>SO<sub>4</sub>.  
Avoid contact with the stop solution. It may cause skin irritations and burns.
  10. **Wash Solution** \*, 1 vial, 30 mL (20X concentrated for 600 mL), pH  $6.5 \pm 0.1$   
see „Preparation of Reagents“.
- \* contain non-mercury preservative

#### 4.1.1 Equipment and material required but not provided

- A microtiter plate calibrated reader (450/620 nm  $\pm$  10 nm) (e.g. the DRG Instruments Microtiter Plate Reader)
- Calibrated variable precision micropipettes
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Timer
- Absorbent paper

### 4.2 Storage Conditions and stability of the Kit

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for two months if stored as described above.

### 4.3 Reagent Preparation

Allow all reagents and required number of strips to reach room temperature prior to use.

#### **Wash Solution**

Dilute *Wash Solution 1+19* (e.g. 10 mL + 190 mL) with fresh and germ free redistilled water. This diluted wash solution has a pH value of  $7.2 \pm 0.2$ .

Consumption: ~ 5 mL per determination.

Crystals in the solution disappear by warming up to 37 °C in a water bath. Be sure that the crystals are completely dissolved before use.

*The diluted Wash Solution is stable for 4 weeks at 2 °C to 8 °C.*

### 4.4 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheets.

### 4.5 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

## 5 SPECIMEN COLLECTION AND PREPARATION

Serum or plasma (EDTA-, heparin- or citrate plasma) can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

*Please note:* Samples containing sodium azide should not be used in the assay.

### 5.1 Specimen Collection

#### **Serum:**

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

#### **Plasma:**

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

### 5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 3 days at 2 °C to 8 °C prior to assaying.

Specimens held for a longer time should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

### 5.3 Specimen Dilution

Prior to assaying each patient specimen is first to be diluted with *Sample Diluent*. For the absorption of rheumatoid factor these prediluted samples then have to be incubated with *IgG-RF-Sorbent*

1. Dilute each patient specimen **1+50** with *Sample Diluent*, e.g. 10 µL of specimen + 0.5 mL of *Sample Diluent*. **Mix well.**
2. Mix well the *IgG-RF-Sorbent* before use.
3. Dilute this prediluted sample **1+1** with *IgG-RF-Sorbent* e.g. 60 µL prediluted sample + 60 µL *IgG-RF-Sorbent*. **Mix well**
4. **Let stand at room temperature for at least 15 minutes, up to a maximum of 2 hours and mix well again.**
5. Take 100 µL of these pretreated samples for the ELISA.

**Please note:** Controls are ready for use and must not be diluted!

## 6 ASSAY PROCEDURE

### 6.1 General Remarks

- It is very important to bring all reagents, samples and controls to room temperature before starting the test run!
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- During 37 °C incubation cover microtiter strips with foil to avoid evaporation.

## 6.2 Test Procedure

Prior to commencing the assay, dilute *Wash Solution*, **prepare patient samples as described in point 5.3**, mix well before pipette and establish carefully the **distribution and identification plan** supplied in the kit for all specimens and controls.

1. Select the required number of microtiter strips or wells and insert them into the holder.  
Please allocate at least:
 

1 well (e.g. A1)	for the substrate blank,
1 well (e.g. B1)	for the <i>Neg. Control</i>
3 wells (e.g. C1-E1)	for the <i>Standard 1-3</i>
1 well (e.g. F1)	for the <i>Pos. Control</i> .

 It is left to the user to determine standards, control and patient samples in duplicate.
2. Dispense
 

<b>100 µL</b> of <i>Neg. Control</i>	into well B1
<b>100 µL</b> of <i>Standard 1</i>	into well C1
<b>100 µL</b> of <i>Standard 2</i>	into well D1
<b>100 µL</b> of <i>Standard 3</i>	into well E1
<b>100 µL</b> of <i>Pos. Control</i>	into well F1 and
<b>100 µL</b> of each pre-treated sample	<u>with new disposable tips</u> into appropriate wells.

 Leave well A1 for substrate blank!
3. Cover wells with foil supplied in the kit. Incubate for **60 minutes at 37 °C**.
4. Briskly shake out the contents of the wells.  
Rinse the wells 5 times with diluted *Wash Solution* (**300 µL per well**). Strike the wells sharply on absorbent paper to remove residual droplets.  
**Important note:**  
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
5. Dispense **100 µL Enzyme Conjugate** into each well, **except A1**.
6. Incubate for **30 minutes at room temperature (20 °C to 25 °C)**.  
*Do not expose to direct sun light!*
7. Briskly shake out the contents of the wells.  
Rinse the wells 5 times with diluted *Wash Solution* (300 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.
8. Add **100 µL of Substrate Solution** into all wells.
9. Incubate for **exactly 10 minutes at room temperature (20 °C to 25 °C) in the dark**.
10. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.  
Any blue color developed during the incubation turns into yellow.  
**Note:** Highly positive patient samples can cause dark precipitates of the chromogen!
11. Read the optical density at **450/620 nm** with a microtiter plate reader **within 30 minutes** after adding the *Stop Solution*.

## 6.3 Measurement

**Adjust** the ELISA microplate or microstrip reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for each control and patient sample in the distribution and identification plan.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable **calculate the mean absorbance values** of all duplicates.

## 7 CALCULATION OF RESULTS

### 7.1 Validation of the Test Run

The test run may be considered valid provided the following criteria are met:

<b>Substrate blank in A1:</b>	Absorbance value <b>lower than 0.100</b>
<b>Neg. Control in B1:</b>	Absorbance value <b>lower than 0.200</b>
<b>Standard 1 (Cut-off) in C1 :</b>	Absorbance value <b>between 0.350 – 0.900</b>
<b>Standard 2 in D1 :</b>	Absorbance value <b>between 0.800 – 1.500</b>
<b>Standard 3 in E1 :</b>	Absorbance value <b>between 1.100 – 2.000</b>
<b>Pos. Control in F1 :</b>	Absorbance value <b>between 0.650 – 3.000</b>

### 7.2 Calculation of quantitative Results

In order to obtain **quantitative results in DU/mL** (DU = DRG Units) plot the (mean) absorbance values of the *Neg. Control* and the 3 Standards 1, 2 and 3 on (linear/linear) graph paper in a system of coordinates against their corresponding concentrations (**0, 50, 200 and 400 DU/mL**) and draw a standard calibration curve (absorbance values on the vertical y-axis, concentrations on the horizontal x-axis).

Read results from this standard curve employing the (mean) absorbance values of each patient specimen and control. All suitable computer programs available can be used for automated result reading and calculation. The following mathematical functions can be used: 4 PL (4 Parameter Logistics) curve fit, Linear regression or Point to Point calculation of the standard curve. We use DRG regression program for windows (4 parameter Rodbart regression). If other regression software is used, the obtained values have to be validated by the user.

*NOTE: Values of additionally (1:10, in total 1:1000) diluted patient samples must be multiplied by the appropriate dilution factor in order to obtain correct results! (Dilution: 1:10 = Dilution factor: 10). (See chapter "5.3 Specimen Dilution").*

### 7.3 Interpretation of quantitative Results

Normal value ranges for this ELISA should be established by each laboratory based on its own patient populations in the geographical areas serviced.

The following values should be considered as a guideline:

<b>NEGATIVE</b>	<b>&lt; 45 DU/mL</b>
<b>CUT-OFF VALUE:</b>	<b>50 DU/mL</b>
<b>GREY ZONE (equivocal):</b>	<b>45 - 55 DU/mL</b>
<b>POSITIVE:</b>	<b>&gt; 55 DU/mL</b>

### 7.4 Calculation of qualitative Results

Absorbance value of **Standard 1 (Cut-off) = CO**

*Example: 0.4 = CO*

### 7.5 Interpretation of qualitative Results

<b>NEGATIVE</b>	<b>Mean OD patient &lt; OD CO –10%</b>
<b>GREY ZONE</b>	<b>OD CO -10% ≤ Mean OD patient ≤ OD CO +10%</b> Repeat test 2 - 4 weeks later - with new patient samples Results in the second test again in the grey zone → NEGATIVE
<b>POSITIVE</b>	<b>Mean OD patient &gt; OD CO +10 %</b>

## 8 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Assay Dynamic Range

The range of the assay is between 0.83 - 400 DU/mL.

### 9.2 Specificity of Antigen (Cross Reactivity)

No cross reactivity was found for Herpes-simplex Virus 1 and 2, Varicella zoster Virus and Epstein-Barr Virus (VCA).

### 9.3 Analytical Sensitivity

The analytical sensitivity of the DRG ELISA was calculated by adding 2 standard deviations from the mean of 20 replicate analyses of the negative control and was found to be 0.83 DU/mL (OD<sub>450nm</sub> 0.055).

### 9.4 Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. (Detected by method comparison with Virion/Serion ELISA with three lots of DRG ELISA, 85 samples, therefrom 69 negative samples are assayed.)

It is 100% (for all three DRG production lots).

### 9.5 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. (Detected by method comparison with Virion/Serion ELISA with three lots of DRG ELISA, 85 samples, therefrom 16 positive samples are assayed.)

It is 100% (for all three DRG production lots).

### 9.6 Method Comparison

The DRG ELISA was compared with the Virion/Serion CMV IgM ELISA. 85 serum samples are assayed.

n = 85		Diamed Eurogen ELISA	
		pos.	neg.
DRG ELISA Lot 1	pos.	16	0
	neg.	0	69

Agreement: 100%

**9.7 Reproducibility**

**9.7.1 Intra-assay**

The intra-assay (within-run) precision of the DRG CMV IgM ELISA was determined by 20 x measurements of 12 serum samples covering the whole measuring range.

Sample	Mean Conc. (DU/mL)	Intra-Assay CV (%)	n
1	21.58	9.55	20
2	30.09	9.54	20
3	5.83	9.41	20
4	105.43	6.52	20
5	112.60	6.20	20
6	87.78	9.42	20
7	350.03	3.09	20
8	195.21	5.71	20
9	212.83	5.51	20
10	464.41	3.74	20
11	327.19	4.87	20
12	456.94	5.86	20

**9.7.2 Inter-assay**

The inter-assay variation of the DRG CMV IgM ELISA was determined with 3 samples with 2 production kits in 10 independent runs with 2 replicates per run.

Sample	Mean Conc. (DU/mL)	Inter-Assay CV (%)	n
1	53.41	14.58	40
2	171.41	10.78	40
3	580.32	11.87	40

**9.8 Recovery**

Samples have been spiked by adding 3 solutions with known concentrations in a 1:1 ratio.

The % recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100 (expected value = (endogenous value + added value) / 2; because of a 1:2 dilution of serum with spike material).

	Sample 1	Sample 2	Sample 3
<b>Concentration [DU/mL]</b>	45.21	27.41	8.60
<b>Average Recovery</b>	93.7	94.0	98.3
<b>Range of Recovery [%]</b>	from	88.9	93.1
	to	97.5	95.6
		83.2	110.9

**9.9 Linearity**

Three samples (serum) containing different amounts of analyte were serially diluted with sample diluent and assayed with the DRG ELISA. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

	Sample 1	Sample 2	Sample 3
<b>Concentration [DU/mL]</b>	536.81	454.08	357.45
<b>Average Recovery</b>	93.0	98.8	102.4
<b>Range of Recovery [%]</b>	from	85.2	93.1
	to	102.3	103.0
		94.8	109.6

## 10 LIMITATIONS OF USE

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values. In immunocompromised patients and newborns serological data only have restricted value.

### 10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

## 11 LEGAL ASPECTS

### 11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

### 11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

## 12 REFERENCES / LITERATURE

1. Lentz, E.B. et al., Detection of Antibody to CMV Induced Early Antigens and Comparison with Four Serologic Assays and Presence of Viruria in Blood Donors, *J. of Clin. Microbiol.*, 26, 133-35 (1988)
2. Von Loon, H.M. et al., Direct Enzyme-Linked Immunosorbent Assay that uses Peroxidase-Labelled Antigens for Determination of Immunoglobulin M Antibody to CMV, *Journal of Clinical Microbiology* 13, 416-422 (1981)
3. Matas, A.J., Simmons, R.L., Fryd, D. and Jajarian, J.S., Persistent, Recurrent and Late CMV Infection, *Transplantation Proceedings*, Vol. XIII, Nr. 1, 114-124 (1981)

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

## SHORT INSTRUCTIONS FOR USE

	<p>All reagents and specimens must be allowed to come to room temperature (18-25°C) before use.</p>
	<p>Leave well A1 for substrate Blank. Dispense 100 µL of Standards and Control into appropriate wells.</p>
	<p>Dispense 100 µL of sample into selected wells. <b>(Please note special sample treatment, point 5.3!)</b></p>
	<p>Cover wells with foil. Incubate for <b>60 minutes</b> at 37 °C.</p>
	<p>Briskly shake out the contents of the wells.</p>
	<p>Rinse the wells 5 times with diluted Wash Solution (300 µL per well).</p>
	<p>Strike the wells sharply on absorbent paper to remove residual droplets.</p>
	<p>Dispense 100 µL of Enzyme-Conjugate into each well.</p>
	<p>Incubate for <b>30 minutes</b> at room temperature.</p>
	<p>Briskly shake out the contents of the wells.</p>
	<p>Rinse the wells 5 times with diluted Wash Solution (300 µL per well).</p>
	<p>Strike the wells sharply on absorbent paper to remove residual droplets.</p>
	<p>Add 100 µL of Substrate Solution to each well.</p>
	<p>Incubate for <b>10 minutes</b> at room temperature.</p>
	<p>Stop the reaction by adding 100 µL of Stop Solution to each well.</p>
	<p>Determine the absorbance of each well at 450 nm.</p>