DCM003-14

Ed. 02/2024

# ESTRADIOL ELISA

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## 18. BIBLIOGRAFIA

1. Smy L, Straseski JA. Measuring estrogens in women, men, and children: Recent advances 2012-2017. *Clin Biochem*. 2018 Dec;62:11-23.
2. Rosner W, Hankinson SE, Sluss PM, Vesper HW, Wierman ME. Challenges to the measurement of estradiol: an endocrine society position statement. *J Clin Endocrinol Metab*. 2013 Apr;98(4):1376-87.
3. van Hellemond IEG, Vriens IJH, Peer PGM, Swinkels ACP, *et al.*; Dutch Breast Cancer Research Group (BOOG). Ovarian Function Recovery During Anastrozole in Breast Cancer Patients With Chemotherapy-Induced Ovarian Function Failure. *J Natl Cancer Inst*. 2017 Dec 1;109(12).
4. Katulski K, Slawek S, Czyzyk A, Podfigurna-Stopa A, *et al.* Bone mineral density in women with polycystic ovary syndrome. *J Endocrinol Invest*. 2014 Dec;37(12):1219-24.
5. Cameron DA, Douglas S, Brown JE, Anderson RA. Bone mineral density loss during adjuvant chemotherapy in pre-menopausal women with early breast cancer: is it dependent on oestrogen deficiency? *Breast Cancer Res Treat*. 2010 Oct;123(3):805-14.
6. Decaroli MC, Rochira V. Aging and sex hormones in males. *Virulence*. 2017 Jul 4;8(5):545-570.
7. Basic QC Practices On-line Course; <http://www.Westgard.com>.
8. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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# ESTRADIOL ELISA

for routine analysis

Immunoenzymatic determination of 17 $\beta$ -Estradiol in human serum or plasma

IVD

LOT

See external label

2°C  8°C  $\Sigma = 96$  tests

REF DKO003

## 1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use  
For Laboratory Professional Use

Estradiol ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of 17 $\beta$ -Estradiol in human serum or plasma from an adult population.

## 2. CLINICAL SIGNIFICANCE

Estradiol (also known as 17 $\beta$ -estradiol, oestradiol or E2) is one of the three main estrogens produced in the human body. The other two are Estrone (E1) and Estriol (E3) which are mainly produced in post-menopausal women and during pregnancy respectively. Estradiol is the predominant estrogen synthesised in pre-menopausal women and has a key role in the regulation of fertility in women.

Estradiol is synthesised in the adrenal glands, the ovaries (women) and testes (men). It is produced by the conversion of testosterone to estradiol by the action of the enzyme aromatase. The enzyme aromatase is also responsible for production of estrone by the conversion of androstenedione<sup>1</sup>.

Estrogens regulate the development of secondary sex characteristics; estradiol is also involved in other physiological functions in different organs and systems such as skin, blood vessels, bone, muscle, gastrointestinal tract, brain, lung and pancreas<sup>2</sup> etc. Quantification of estrogens (and androgens) are useful in the assessment and management of different sex-hormone related disorders including hypogonadism, hirsutism, amenorrhea, fertility and ovarian tumours.

Due to the importance of estradiol in the regulation of the menstrual cycle, low levels are associated with disturbance of the cycle and amenorrhea. Measurements can therefore be used to assess ovarian functionality and support diagnosis of ovarian function recovery in women treated with aromatase inhibitors<sup>3</sup>.

Estradiol levels are known to be low in men where its biological function is not clearly understood. There is some indication that it is involved in several physiological and pathological conditions such as prostate cancer and

infertility. Conversely, elevated levels in men, due to high activity of aromatase, inhibits the synthesis of luteinising hormone from the pituitary gland leading to a decrease in testosterone and onset of hypogonadism and gynecomastia.

Estradiol has a key role in bone homeostasis where low levels are associated with loss of bone mass. It has influence on osteoblast and osteoclast activity to support the maintenance of bone homeostasis and is considered a good predictor of bone mass<sup>1,2,4-6</sup>.

## 3. PRINCIPLE OF THE METHOD

The Estradiol ELISA is a competitive enzyme immunometric assay (ELISA) where 17 $\beta$ -Estradiol (antigen) in the sample competes with the antigenic 17 $\beta$ -Estradiol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti 17 $\beta$ -Estradiol coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added. The colour intensity is inversely proportional to the 17 $\beta$ -Estradiol concentration in the sample.

17 $\beta$ -Estradiol concentration in the sample is calculated through a calibration curve.

## 4. REAGENTS, MATERIALS AND INSTRUMENTATION

### 4.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 0.5 mL each)

Human serum, ProClin >0.0015%, NaN<sub>3</sub> <0.1% and BSA

CAL0	REF DCE002/0306-0
CAL1	REF DCE002/0307-0
CAL2	REF DCE002/0308-0
CAL3	REF DCE002/0309-0
CAL4	REF DCE002/0310-0
CAL5	REF DCE002/0311-0

2. Control (1 vial, 0.5 mL)

Human serum, ProClin >0.0015%, NaN<sub>3</sub> <0.1% and BSA

Concentration of Control is indicated on the Certificate of Analysis

REF DCE045/0303-0

3. Conjugate (1 vial, 22 mL)  
17β-Estradiol conjugated with horseradish peroxidase (HRP). ProClin <0.0015% and BSA  
**REF DCE002/0302-0**
4. Coated Microplate (1 microplate breakable)  
Anti 17β-Estradiol antibodies adsorbed on microplate  
**REF DCE002/0303-0**
5. TMB Substrate (1 vial, 15 mL)  
H<sub>2</sub>O<sub>2</sub>-TMB (0.26 g/L) (avoid any skin contact)  
ProClin <0.0015%  
**REF DCE004-0**
6. Stop Solution (1 vial, 15 mL)  
Sulphuric acid 0.15M (avoid any skin contact)  
**REF DCE005-0**
7. 10X Conc. Wash Solution (1 vial, 50 mL)  
Phosphate buffer 0.2M pH 7.4. Contains ProClin >0.0015%  
**REF DCE054-0**

#### 4.2. Materials required but not provided

Distilled water

#### 4.3. Auxiliary materials and instrumentation

Automatic dispenser  
Precision Pipetting Devices  
Microplate reader (450 nm, 620-630 nm)

#### 5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, control and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]  
Skin sensitivity, Category 1

Warning

Contains: ProClin 300

#### Hazard statements:

H317 - May cause an allergic skin reaction.

#### Precautionary statements:

P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.  
P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.  
P321 - Specific treatment (see supplemental first aid instruction on this label).  
P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.  
P362+P364 - Take off contaminated clothing and wash it before reuse.

- Some reagents (calibrators and control) contain small amounts of Sodium Azide (NaN<sub>3</sub>) <0.1% as preservative. Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through large with amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to direct sunlight, metals or oxidants. Do not freeze the solution.

#### 6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

## 7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

## 8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	24 hours
Freeze/thaw cycles	1 cycle

## 9. PROCEDURE

### 9.1. Preparation of Calibrators and Controls

The Calibrators are ready to use and have the following concentration of 17β-Estradiol:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
pg/mL	0	20	120	300	600	2000

The calibrators are ready to use; the concentration of the Calibrator is printed on the label.

The control is ready to use.

### 9.2. Preparation of the Conjugate

The conjugate is ready to use. Mix gently for 5 minutes with a rotating mixer.

### 9.3. Preparation of the Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals

completely by rinsing of the bottle, then mix until crystals are completely dissolved.

### 9.4. Preparation of Samples

The determination of 17β-Estradiol can be performed in human serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Fasting samples are not necessary and no special sample preparations are required.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation. Testing of heat-inactivated sera is not recommended.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Avoid repetitive freezing and thawing of samples. Before using, mix gently for 5 minutes with a rotating mixer.

### 9.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	25 µL		
Sample/Control		25 µL	
Conjugate	200 µL	200 µL	

Incubate 2 h at +37°C.

Remove the contents from each well, wash the wells three times with 300 µL of diluted wash solution.

**Important note:** during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

**Automatic washer:** if you use automated equipment, wash the wells at least 5 times.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at room temperature (22÷28°C) for 30 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

---

Shake gently the microplate.  
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

---

## 10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories<sup>7</sup>.

## 11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is recommended.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

## Conversion of units

To convert results to SI units:  
 $\text{pmol/L} = \text{pg/mL} \times 3.671$

## 12. MEASURING RANGE

The assay measuring range (AMR) is 29.6 – 2000 pg/mL (108.7 – 7342 pmol/L).

Any value that reads below 29.6 pg/mL (108.7 pmol/L) should be reported as “< 29.6 pg/mL (108.7 pmol/L)”. Any value that reads above 2000 pg/mL (7342 pmol/L) should be reported as “> 2000 pg/mL (7342 pmol/L)”.

## 13. METROLOGY AND TRACEABILITY

The calibrators of this kit are traceable to the Estradiol standard from the National Metrology Institute of Japan (6004-a).

## 14. EXPECTED VALUES

The following ranges were determined using the Estradiol ELISA and are provided for information only. The 90 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

Adults	n	Median (pg/mL)	Reference interval (pg/mL)
<b>Males</b>	120	34.3	<29.6 – 64.2
<b>Females</b>			
Luteal phase	40	71.6	<29.6 – 212.4
Follicular phase	39	46.5	<29.6 – 135.6
Ovulatory phase	41	57.1	<29.6 – 273.5
Post-menopausal	120	32.9	<29.6 – 65.8
<b>Pregnancy</b>			
1 <sup>st</sup> trimester	40	1120.1	301.2 – >2000
2 <sup>nd</sup> trimester	39	>2000	1850.3 – >2000
3 <sup>rd</sup> trimester	40	>2000	>2000

Children	n	Median (pg/mL)	Reference interval (pg/mL)
<b>Females</b>			
1 – 10 years	39	56.7	<29.6 – 115.1
10 – 15 years	41	79.8	46.9 – 155.5
15 – 19 years	40	138.3	88.8 – 264.1
<b>Males</b>			
1 – 10 years	40	58.5	<29.6 – 115.7
10 – 15 years	39	39.2	<29.6 – 105.1
15 – 19 years	39	31.7	<29.6 – 91.3

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

## 15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

### 15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	6.8 pg/mL
Limit of Detection (LoD)	14.6 pg/mL
Limit of Quantitation (LoQ)	29.6 pg/mL

### 15.2. Trueness

Trueness has been demonstrated through a recovery test for the Estradiol ELISA with the Estradiol standard from the National Metrology Institute of Japan (code 6004-a).

### 15.3. Precision

Precision of the Estradiol ELISA was determined by performing a complex precision study.

**Repeatability:** A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (pg/mL)	Within run (Repeatability)	
			SD	CV%
1	75	51.7	6.1	11.9%
2	75	134.2	11.3	8.4%
3	75	264.0	25.4	9.6%
4	75	627.7	56.4	9.0%
5	75	838.1	61.5	7.3%
6	75	1346.6	120.2	8.9%

**Reproducibility:** A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Results for the combined data from 2 lots is shown below:

Sample	n	Mean Conc. (pg/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	52.4	8.2	15.7%
2	150	139.0	21.2	15.3%
3	150	274.7	39.0	14.2%
4	150	633.9	66.0	10.4%
5	150	860.0	95.3	11.1%
6	150	1374.7	159.8	11.6%

### 15.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For estradiol concentration by Estradiol ELISA, the measurement procedure shows linearity for the interval from 16.5 to 2271.3pg/mL (60.57 to 8337.9 pmol/L) within the allowable deviation of linearity (ADL) of  $\pm 15\%$ .

### 15.5. Method comparison

The Estradiol ELISA was compared against a commercially available quantitative manual assay, following CLSI EP-09C, "Method Comparison and Bias Estimation Using Patient Samples". A total of 106 samples, selected to represent a wide range of estradiol concentrations, was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
106	0.98 [0.93 to 1.04]	14.5 [8.6 to 18.2]	0.99

### 15.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross reactant	Concentration tested	Mean % Cross reactivity
Estradiol	100 ng/mL	1.0%
Estrone	100 ng/mL	0.9%
Fulvestrant	100 ng/mL	0.4%
Testosterone	100 ng/mL	0.0%
Cortisol	1000 ng/mL	0.0%
Progesterone	200 ng/mL	0.0%
DHEA-S	10000 ng/mL	0.0%
17OH Progesterone	100 ng/mL	0.0%
5 $\alpha$ -dihydrotestosterone (DHT)	500 ng/mL	0.0%
Androstenedione	100 ng/mL	0.0%
Androsterone	100 ng/mL	0.0%
Cortisone	1000 ng/mL	0.0%
Prednisolone	200 ng/mL	0.0%
Danazol	1000 ng/mL	0.0%
2-Methoxyestradiol	10 ng/mL	0.7%
17 $\alpha$ -Ethinylestradiol	200 ng/mL	0.0%
17 $\beta$ -Estradiol 17-sulphate	100 ng/mL	0.0%

The following substances do not interfere with a bias of > ±15% in the Estradiol ELISA assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	3.89 mg/dL
Bilirubin, unconjugated	20 mg/dL
Haemoglobin	640 mg/dL
Total Protein	8.3 g/dL
Triglycerides	2000 mg/dL

### 15.7. Serum-plasma study

The Estradiol ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP-35 guidelines. A total of 22 samples (18 native, 4 spiked) to cover the range were evaluated. Passing-Bablok regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
SST	1.03 [0.98 to 1.07]	-3.5 [-37.2 to 30.2]	1.00
Lithium Heparin	1.04 [0.97 to 1.12]	-0.8 [-59.0 to 57.4]	0.99
Sodium Heparin	1.10 [1.03 to 1.17]	-8.05 [-62.2 to 46.1]	0.99
EDTA	1.06 [1.01 to 1.11]	-7.6 [-46.3 to 31.1]	0.99

### 16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays<sup>8</sup>. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

### 17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

### 18. BIBLIOGRAPHY

1. Smy L, Straseski JA. Measuring estrogens in women, men, and children: Recent advances 2012-2017. *Clin Biochem.* 2018 Dec;62:11-23.
2. Rosner W, Hankinson SE, Sluss PM, Vesper HW, Wierman ME. Challenges to the measurement of estradiol: an endocrine society position statement. *J Clin Endocrinol Metab.* 2013 Apr;98(4):1376-87.
3. van Hellemond IEG, Vriens IJH, Peer PGM, Swinkels ACP, *et al.*; Dutch Breast Cancer Research Group (BOOG). Ovarian Function Recovery During Anastrozole in Breast Cancer Patients With Chemotherapy-Induced Ovarian Function Failure. *J Natl Cancer Inst.* 2017 Dec 1;109(12).
4. Katulski K, Slawek S, Czyzyk A, Podfigurna-Stopa A, *et al.* Bone mineral density in women with polycystic ovary syndrome. *J Endocrinol Invest.* 2014 Dec;37(12):1219-24.
5. Cameron DA, Douglas S, Brown JE, Anderson RA. Bone mineral density loss during adjuvant chemotherapy in pre-menopausal women with early breast cancer: is it dependent on oestrogen deficiency? *Breast Cancer Res Treat.* 2010 Oct;123(3):805-14.
6. Decaroli MC, Rochira V. Aging and sex hormones in males. *Virulence.* 2017 Jul 4;8(5):545-570.
7. Basic QC Practices On-line Course; <http://www.Westgard.com>.
8. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

### 19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

### 20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: [www.diametra.com](http://www.diametra.com).

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1. Smy L, Straseski JA. Measuring estrogens in women, men, and children: Recent advances 2012-2017. *Clin Biochem*. 2018 Dec;62:11-23.
2. Rosner W, Hankinson SE, Sluss PM, Vesper HW, Wierman ME. Challenges to the measurement of estradiol: an endocrine society position statement. *J Clin Endocrinol Metab*. 2013 Apr;98(4):1376-87.
3. van Hellemond IEG, Vriens IJH, Peer PGM, Swinkels ACP, *et al.*; Dutch Breast Cancer Research Group (BOOG). Ovarian Function Recovery During Anastrozole in Breast Cancer Patients With Chemotherapy-Induced Ovarian Function Failure. *J Natl Cancer Inst*. 2017 Dec 1;109(12).
4. Katulski K, Slawek S, Czyzyk A, Podfigurna-Stopa A, *et al.* Bone mineral density in women with polycystic ovary syndrome. *J Endocrinol Invest*. 2014 Dec;37(12):1219-24.
5. Cameron DA, Douglas S, Brown JE, Anderson RA. Bone mineral density loss during adjuvant chemotherapy in pre-menopausal women with early breast cancer: is it dependent on oestrogen deficiency? *Breast Cancer Res Treat*. 2010 Oct;123(3):805-14.
6. Decaroli MC, Rochira V. Aging and sex hormones in males. *Virulence*. 2017 Jul 4;8(5):545-570.

7. Basic QC Practices On-line Course; <http://www.Westgard.com>.
8. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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