Instructions for Use

Filariasis ELISA

IVD

REF EIA-5972

Σ 96

DRG Instruments GmbH, Germany

Distributed by:

DRG International, Inc., USA

**Please use only the valid version of the Instructions for Use provided with the kit.
 Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
 Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
 Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
 Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.
 Por favor, usar a versão válida das instruções de utilização fornecidas com o kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

1	INTRODUCTION.....	3
2	INTENDED USE.....	3
3	PRINCIPLE OF THE ASSAY	3
4	MATERIALS.....	4
5	STABILITY AND STORAGE	4
6	REAGENT PREPARATION	5
7	SAMPLE COLLECTION AND PREPARATION.....	5
8	ASSAY PROCEDURE	6
9	RESULTS.....	7
10	SPECIFIC PERFORMANCE CHARACTERISTICS	7
11	LIMITATIONS OF THE PROCEDURE.....	8
12	PRECAUTIONS AND WARNINGS	8
1	EINLEITUNG.....	9
2	VERWENDUNGSZWECK.....	9
3	TESTPRINZIP	9
4	MATERIALIEN	10
5	STABILITÄT UND LAGERUNG	10
6	VORBEREITUNG DER REAGENZIEN.....	11
7	ENTNAHME UND VORBEREITUNG DER PROBEN	11
8	TESTDURCHFÜHRUNG	12
9	BERECHNUNG DER ERGEBNISSE	13
10	TESTMERKMALE	13
11	GRENZEN DES VERFAHRENS.....	14
12	SICHERHEITSMASSNAHMEN UND WARNHINWEISE	14
1	INTRODUZIONE.....	15
2	USO PREVISTO	15
3	PRINCIPIO DEL TEST	15
4	MATERIALI	16
5	MODALITÀ DI CONSERVAZIONE	16
6	PREPARAZIONE DEI REAGENTI.....	17
7	PRELIEVO E PREPARAZIONE DEI CAMPIONI	17
8	PROCEDIMENTO	18
9	RISULTATI.....	19
10	CARATTERISTICHE DEL TEST.....	19
11	LIMITAZIONI	20
12	PRECAUZIONI E AVVERTENZE.....	20

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1 INTRODUCTION

Filariasis (or philariasis) is a parasitic disease caused by an infection with roundworms of the Filarioidea type. These are spread by blood-feeding black flies and mosquitoes. This disease belongs to the group of diseases called helminthiases. Eight known filarial nematodes use humans as their definitive hosts. These are divided into three groups according to the niche within the body they occupy:

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Wuchereria bancrofti</i> , <i>Brugia malayi</i> , and <i>Brugia timori</i>	Lymphatic filariasis	Recurring lymphadenitis, Fever, Disfigurement of the limbs (elephantiasis) and genitalia (hydroceles and other anatomical changes in the male genitalia) by worms occupying the lymphatic system	Mosquito vectors include <i>Mansonia</i> , <i>Anopheles</i> , and <i>Aedes</i> mosquitoes
<i>Loa loa</i> (the eye worm), <i>Mansonella streptocerca</i> , and <i>Onchocerca volvulus</i>	Subcutaneous filariasis (Loiasis; Streptocerciasis; Onchocerciasis)	Calabar swellings (transient localized angioedema) and subconjunctival migration of the worm (eye worm); chronic, pruritic dermatitis and lymphadenopathy; dermatitis, keratitis, chorioretinitis, blindness (river blindness)	Female black fly of the genus <i>Simulium</i>
<i>Mansonella perstans</i> and <i>Mansonella ozzardi</i>	Serous cavity filariasis (Perstans filariasis)	transient angioedema, pruritus of arms, face and other parts of the body, fever, headache, arthralgias; <i>M. ozzardi</i> infections are usually asymptomatic	Blood-sucking flies called midges (<i>Culicoides grahamsi</i> and <i>C. austeni</i>)

The adult worms, which usually stay in one tissue, release early larval forms known as microfilariae into the host's bloodstream. These circulating microfilariae can be taken up with a blood meal by the arthropod vector; in the vector, they develop into infective larvae that can be transmitted to a new host.

Infection or presence of pathogen may be identified by:

- Microscopy
- PCR
- Serology: e.g. ELISA

2 INTENDED USE

The Filariasis ELISA is intended for the qualitative determination of specific antibodies against filariae in human serum or plasma (citrate or heparin).

3 PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4 MATERIALS

4.1 Reagents supplied

- **Microtiterplate:**
12 break-apart 8-well snap-off strips coated with filarial antigens; in resealable aluminium foil.
- **Sample Dilution Buffer:**
1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:**
1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):**
1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M); pH 7.2 ± 0.2; for washing the wells; white cap.
- **Conjugate:**
1 bottle containing 20 ml of peroxidase labelled Protein A/G; coloured blue, ready to use; black cap; ≤ 0.02% (v/v) MIT.
- **TMB Substrate Solution:**
1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:**
1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:**
1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:**
1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2 Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3 Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5 STABILITY AND STORAGE

Store the kit at 2 °C - 8 °C.

The opened reagents are stable up to the expiry date stated on the label when stored at 2 °C - 8 °C.

6 REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20 °C - 25 °C) and mix them before starting the test run!

6.1 Microtiterplate

The break-apart snap-off strips are coated with filarial antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2 °C - 8 °C.

6.2 Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water.

The diluted buffer is stable for 5 days at room temperature (20 °C - 25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

6.3 TMB Substrate Solution

The reagent is ready to use and has to be stored at 2 °C - 8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7 SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay.

If the assay is performed within 5 days after sample collection, the samples should be kept at 2 °C - 8 °C; otherwise they should be aliquoted and stored deep-frozen (-70 °C to -20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

7.1 Sample Dilution

Before assaying, all samples should be diluted **1+100** with Sample Dilution Buffer.

Dispense 10 µL sample and 1 mL Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8 ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 °C ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 °C ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 seconds. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20 °C - 25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1 Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9 RESULTS

9.1 Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < **0.100**
- **Negative Control:** Absorbance value < **0.200** and < **Cut-off**
- **Cut-off Control:** Absorbance value **0.150 – 1.300**
- **Positive Control:** Absorbance value > **Cut-off**

If these criteria are not met, the test is not valid and must be repeated.

9.2 Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86/2 = 0.43
 Cut-off = 0.43

9.2.1 Results in Units [DU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{DRG Units} = \text{DU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ DU (Units)}$

9.3 Interpretation of Results

Cut-off	10 DU	-
Positive	> 11 DU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 DU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 DU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.

10 SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact DRG.

10.1 Precision

<u>Intra-assay</u>	<u>n</u>	<u>Mean (E)</u>	<u>CV (%)</u>
#1	12	0.435	2.07
#2	12	0.334	3.13
#3	12	0.143	3.97

<u>Inter-assay</u>	<u>n</u>	<u>Mean (DU)</u>	<u>CV (%)</u>
#1	10	11.85	1.81
#2	10	9.45	2.67
#3	10	11.91	2.60

10.2 Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 99.1% (95% confidence interval: 94.99% - 99.98%).

10.3 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 94.7% (95% confidence interval: 82.25% - 99.36%).

10.4 Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5 Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters revealed only minimal cross-reactivity to *Fasciola* specimens.

However, in endemic areas, double infection as well as past infection with other parasites and worms should be considered.

11 LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12 PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1 Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet.

12.2 Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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13 BIBLIOGRAPHY / LITERATUR / BIBLIOGRAFIA / BIBLIOGRAFÍA

1. Alli, R.; Kulkarni, S.; Reddy, M. V.; Harinath, B. C. (2001): Evaluation of sevafilachek immunoassays and rapid ICT-filariasis test for detection of bancroftian filariasis. In *Indian journal of clinical biochemistry : IJCB* 16 (2), pp. 207–210. DOI: 10.1007/BF02864863.
2. Burbelo, Peter D.; Ramanathan, Roshan; Klion, Amy D.; Iadarola, Michael J.; Nutman, Thomas B. (2008): Rapid, novel, specific, high-throughput assay for diagnosis of Loa loa infection. In *Journal of Clinical Microbiology* 46 (7), pp. 2298–2304. DOI: 10.1128/JCM.00490-08.
3. Chen, Jun-Hu; Wang, Hen; Chen, Jia-Xu; Bergquist, Robert; Tanner, Marcel; Utzinger, Jurg; Zhou, Xiao-Nong (2012): Frontiers of parasitology research in the People's Republic of China: infection, diagnosis, protection and surveillance. In *Parasites & vectors* 5, p. 221. DOI: 10.1186/1756-3305-5-221.
4. Cuello, M. Rivera; Cuadros, E. Nunez; Claros, A. Medina; Hortelano, M. Garcia; Fontelos, P. Martin; Pena, M. J. Mellado (2009): Filariasis en pacientes procedentes de area endemica. Presentacion de una serie de 14 casos. In *Anales de pediatria (Barcelona, Spain : 2003)* 71 (3), pp. 189–195. DOI: 10.1016/j.anpedi.2009.04.022.
5. Freedman, David O. (2006): Onchocerciasis. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): *Tropical infectious diseases. Principles, pathogens & practice*. 2nd ed. Philadelphia: Churchill Livingstone, pp. 1176–1188.
6. Hotez, Peter J.; Bottazzi, Maria Elena; Strych, Ulrich; Chang, Li-Yen; Lim, Yvonne A. L.; Goodenow, Maureen M.; AbuBakar, Sazaly (2015): Neglected tropical diseases among the Association of Southeast Asian Nations (ASEAN): overview and update. In *PLoS neglected tropical diseases* 9 (4), e0003575. DOI: 10.1371/journal.pntd.0003575.
7. Joseph, Hayley M.; Melrose, Wayne (2010): Applicability of the Filter Paper Technique for Detection of Antifilarial IgG(4) Antibodies Using the Bm14 Filariasis CELISA. In *Journal of parasitology research* 2010. DOI: 10.1155/2010/594687.
8. Klion, Amy D.; Nutman, Thomas B. (2006): Loiasis and Mansonella Infections. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): *Tropical infectious diseases. Principles, pathogens & practice*. 2nd ed. Philadelphia: Churchill Livingstone, pp. 1163–1175.
9. Kyelem, Dominique; Fischer, Peter U.; Brattig, Norbert W. (2011): The diagnostics and control of neglected tropical helminth diseases. Foreword. In *Acta tropica* 120 Suppl 1, S1-3. DOI: 10.1016/j.actatropica.2011.04.001.
10. Lloyd, Melanie M.; Gilbert, Rebecca; Taha, Nathalie Tebao; Weil, Gary J.; Meite, Aboulaye; Kouakou, Ilunga M. M.; Fischer, Peter U. (2015): Conventional parasitology and DNA-based diagnostic methods for onchocerciasis elimination programmes. In *Acta tropica* 146, pp. 114–118. DOI: 10.1016/j.actatropica.2015.03.019.
11. Nutman, Thomas B.; Kazura, James W. (2006): Filariasis. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): *Tropical infectious diseases. Principles, pathogens & practice*. 2nd ed. Philadelphia: Churchill Livingstone, pp. 1152–1162.
12. Rebollo, Maria P.; Bockarie, Moses J. (2013): Rapid diagnostics for the endgame in lymphatic filariasis elimination. In *The American journal of tropical medicine and hygiene* 89 (1), pp. 3–4. DOI: 10.4269/ajtmh.13-0202.
13. Steel, Cathy; Golden, Allison; Stevens, Eric; Yokobe, Lindsay; Domingo, Gonzalo J.; los Santos, Tala de; Nutman, Thomas B. (2015): Rapid Point-of-Contact Tool for Mapping and Integrated Surveillance of Wuchereria bancrofti and Onchocerca volvulus Infection. In *Clinical and vaccine immunology : CVI* 22 (8), pp. 896–901. DOI: 10.1128/CVI.00227-15.
14. Tambo, Ernest; Ai, Lin; Zhou, Xia; Chen, Jun-Hu; Hu, Wei; Bergquist, Robert et al. (2014): Surveillance-response systems: the key to elimination of tropical diseases. In *Infectious diseases of poverty* 3, p. 17. DOI: 10.1186/2049-9957-3-17.

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14 SCHEME OF THE ASSAY

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all samples and standards/controls on the plate layout supplied in the kit.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Negative Control	Cut-off Control	Positive Control	Sample (diluted 1+100)
Negative Control	-	100 µL	-	-	-
Cut-off Control	-	-	100 µL	-	-
Positive Control	-	-	-	100 µL	-
Sample (diluted 1+100)	-	-	-	-	100 µL
Cover wells with foil supplied in the kit Incubate for 1 h at 37 °C Wash each well three times with 300 µL of Washing Buffer					
Conjugate	-	100 µL	100 µL	100 µL	100 µL
Incubate for 30 min at room temperature (20 °C - 25 °C) Do not expose to direct sunlight Wash each well three times with 300 µL of Washing Buffer					
TMB Substrate	100 µL	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark					
Stop Solution	100 µL	100 µL	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	Numero di Catalogo	Número de catálogo	Référence de catalogue
	Batch code *	Chargencode *	Codice del lotto	Codigo de lote	Numéro de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservacion	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Caution *	Achtung *			
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité
MTP	Microtiterplate	Mikrotiterplatte	Piastre di Microtitolazione	Placa de Microtitulación	Plaque de Microtitrage
CONJ	Conjugate	Konjugat	Coniugato	Conjugado	Conjugué
CONTROL -	Negative Control	Negativkontrolle	Controllo negativo	Control positivo	Contrôle négatif
CONTROL +	Positive Control	Positivkontrolle	Controllo positivo	Control positivo	Contrôle positif
CUT OFF	Cut off control	Cut off-Kontrolle	Controllo cut-off	Control cut-off	Contrôle cut-off
CAL	Standard or Calibrator	Standard oder Kalibrator	Standard o Calibratore	Estándar o Calibrador	Standard o Etalon
DIL	Sample Dilution Buffer	Probenverdünnungspuffer	Tampone di Diluizione del Campione	Tampón de Dilución de Muestra	Tampon de Dilution d'Échantillon
DIL A	IgA Sample Dilution Buffer	IgA-Probenverdünnungspuffer	Tampone di Diluizione del Campione IgA	Tampón de Dilución de Muestra IgA	Tampon de Dilution d'Échantillon IgA
DIL G	IgG Sample Dilution Buffer	IgG-Probenverdünnungspuffer	Tampone di Diluizione del Campione IgG	Tampón de Dilución de Muestra IgG	Tampon de Dilution d'Échantillon IgG
DIL M	IgM Sample Dilution Buffer	IgM-Probenverdünnungspuffer	Tampone di Diluizione del Campione IgM	Tampón de Dilución de Muestra IgM	Tampon de Dilution d'Échantillon IgM
SOLN STOP	Stop solution	Stopplösung	Soluzione bloccante	Solución de parada	Solution d'arrêt
SUB TMB	TMB Substrate solution	TMB-Substratlösung	Soluzione substrato TMB	Solución de sustrato de TMB	Solution de substrat TMB
WASH BUF 20x	Washing Buffer 20x concentrated	Waschpuffer 20x konzentriert	Tampone di lavaggio concentrazione x20	Tampón de lavado concentrado x20	Tampon de lavage concentré 20 x