DCM006-11

Ed. 03/2024

# PROGESTERONE ELISA

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#### 18. BIBLIOGRAFIA

1. Wu SP, Li R, DeMayo FJ. Progesterone Receptor Regulation of Uterine Adaptation for Pregnancy. *Trends Endocrinol Metab.* 2018 Jul;29(7):481-491.
2. Kasum M, Simunić V, Vrčić H, Stanić P, Orešković S, Beketić-Orešković L. Follicular progesterone elevations with ovulation induction for IVF. *Gynecol Endocrinol.* 2014 Aug;30(8):537-41.
3. Di Renzo GC, Giardina I, Clerici G, Brillo E, Gerli S. Progesterone in normal and pathological pregnancy. *Horm Mol Biol Clin Investig.* 2016 Jul 1;27(1):35-48.
4. Feng C, Chen ZY, Zhang J, Xu H, Zhang XM, Huang XF. Clinical utility of serum reproductive hormones for the early diagnosis of ectopic pregnancy in the first trimester. *J Obstet Gynaecol Res.* 2013 Feb;39(2):528-35.
5. Rausch ME, Sammel MD, Takacs P, Chung K, Shaunik A, Barnhart KT. Development of a multiple marker test for ectopic pregnancy. *Obstet Gynecol.* 2011 Mar;117(3):573-82.
6. Basic QC Practices On-line Course; <http://www.Westgard.com>.
7. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

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# PROGESTERONE ELISA

for routine analysis

Direct immunoenzymatic determination of Progesterone in human serum or plasma

IVD

LOT

See external label

2°C 8°C

 $\Sigma = 96$  tests

REF DKO006

## 1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use  
For Laboratory Professional Use

Progesterone ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of progesterone in human serum or plasma. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis and monitoring of disorders of the ovaries and placenta.

## 2. CLINICAL SIGNIFICANCE

Progesterone is a hormone synthesised in the ovaries, adrenal glands and placenta. It diffuses into circulation where it binds to plasma steroid-binding proteins such as albumin and corticosteroid binding protein. Progesterone is then transported toward the target organs including the hypothalamus, ovaries, uterus but also the liver, where it is catabolised<sup>1</sup>.

Progesterone is fundamental in the early phase of pregnancy since it is required for uterine receptivity for embryo implantation<sup>1</sup>. During the early follicular phase of the female menstrual cycle, the majority of circulating progesterone is of adrenal production. However, in the late follicular phase, there is a rise in progesterone serum levels released by the corpus luteum, a luteinised ovarian follicle<sup>1</sup>. This physiological increase of progesterone is crucial to ensure the secretory function of endometrium and a normal embryo implantation<sup>2</sup>. It does so by controlling the maternal immune system to promote maternal immunotolerance to avoid embryo rejection and reduces uterine contractility and improves utero-placental circulation for proper growth of the embryo<sup>3</sup>.

Measurement of progesterone is useful to aid the assessment of ovulation and to aid in the diagnosis of ovarian dysfunctions that could lead to infertility. In addition, assessment of progesterone is also used as an aid in the diagnosis of ectopic pregnancy (EP), where levels are lower in EP than intrauterine pregnancy<sup>4,5</sup>.

Progesterone levels are relatively low in children and postmenopausal women. Adult males have levels similar to

those in women during the follicular phase of the menstrual cycle.

## 3. PRINCIPLE OF THE METHOD

The Progesterone ELISA is a competitive enzyme immunometric assay (ELISA) where progesterone (antigen) in the sample competes with the antigenic progesterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti progesterone coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate ( $H_2O_2$ ) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution ( $H_2SO_4$ ) is added. The colour intensity is inversely proportional to the progesterone concentration in the sample.

Progesterone concentration in the sample is calculated through a calibration curve.

## 4. REAGENTS, MATERIALS AND INSTRUMENTATION

### 4.1. Reagents and materials supplied in the kit

#### 1. Calibrators (5 vials, 1 mL each)

ProClin &gt;0.0015%

CAL0

REF DCE002/0606-0

CAL1

REF DCE002/0607-0

CAL2

REF DCE002/0608-0

CAL3

REF DCE002/0609-0

CAL4

REF DCE002/0610-0

CAL5

REF DCE002/0611-0

#### 2. Control (1 vial, 1 mL)

ProClin &gt;0.0015%

Concentration of Control is indicated on the Certificate of Analysis

REF DCE045/0603-0

#### 3. Conjugate (1 vial, 22 mL)

Progesterone conjugated with horseradish peroxidase (HRP), ProClin >0.0015% and BSA 0.1%

REF DCE002/0602-0

#### 4. Coated Microplate (1 microplate breakable)

Antibody anti Progesterone coated on the microplate

REF DCE002/0603-0

5. **TMB Substrate** (1 vial, 15 mL)  
 H<sub>2</sub>O<sub>2</sub>-TMB 0.26 g/L (avoid any skin contact)  
 ProClin <0.0015% **REF DCE004-0**
6. **Stop Solution** (1 vial, 15 mL)  
 Sulphuric acid 0.15 mol/L (avoid any skin contact)  
**REF DCE005-0**
7. **10X Conc. Wash Solution** (1 vial, 50 mL)  
 Phosphate buffer 0.2M pH 7.4, ProClin >0.0015%  
**REF DCE054-0**

#### 4.2. Materials required but not provided

Distilled water

#### 4.3. Auxiliary materials and instrumentation

Automatic dispenser  
 Precision Pipetting Devices  
 Microplate reader (450 nm, 620-630 nm)

#### 5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, controls, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to direct sunlight, metals or oxidants. Do not freeze the solution.

#### 6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

#### 7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

#### 8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	96 hours
Freeze/thaw cycles	3 cycles

## 9. PROCEDURE

### 9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer. The calibrators are ready for use and have the following concentration:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
<b>ng/mL</b>	0	0.2	1.0	5.0	15.0	40.0

The control is ready to use; the concentration is printed on the label.

### 9.2. Preparation of the Conjugate

The conjugate is ready to use.

### 9.3. Preparation of the Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

### 9.4. Preparation of Samples

The determination of progesterone can be performed in human or plasma samples.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

### 9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for the Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	<b>20 µL</b>		
Sample/ Control		<b>20 µL</b>	

Incubate at 25°C ±0.5°C for 5 minutes.

Conjugate	<b>200 µL</b>	<b>200 µL</b>	
-----------	---------------	---------------	--

Incubate 1 h at +25°C ±0.5°C.

Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

**Important note:** during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

**Automatic washer:** if you use automated equipment, wash the wells at least 5 times.

TMB Substrate	<b>100 µL</b>	<b>100 µL</b>	<b>100 µL</b>
---------------	---------------	---------------	---------------

Incubate at 22 – 28°C for 15 minutes in the dark.

Stop Solution	<b>100 µL</b>	<b>100 µL</b>	<b>100 µL</b>
---------------	---------------	---------------	---------------

Shake gently the microplate.

Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

## 10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of the control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give

control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories<sup>6</sup>.

## 11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

## 12. MEASURING RANGE

The assay measuring range (AMR) is 0.26 – 40 ng/mL. Any value that reads below 0.26 ng/mL should be reported as “< 0.26 ng/mL”. Any value that reads above 40 ng/mL should be reported as “> 40 ng/mL”.

## 13. METROLOGY AND TRACEABILITY

The calibrators of this kit are traceable to the Progesterone standard from the National Metrology Institute of Japan (6003-a).

## 14. EXPECTED VALUES

The following ranges were determined using the Progesterone ELISA and are provided for information only. The 90 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

Adults	n	Median (ng/mL)	Expected values (ng/mL)
<b>Males</b>	116	0.8	<0.26 - 1.77
<b>Females</b>			
<b>Follicular phase</b>	118	0.434	< 0.26 - 1.94
<b>Luteal phase</b>	42	4.72	2.16 - 18.55
<b>Postmenopausal</b>	120	<0.26	<0.26 - 0.82

<b>Pregnancy</b>			
<b>1<sup>st</sup> trimester</b>	39	37.09	14.80 - >40.00
<b>2<sup>nd</sup> trimester</b>	39	>40.00	38.47 - >40.00
<b>3<sup>rd</sup> trimester</b>	40	>40.00	>40.00

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

## 15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

### 15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.02 ng/mL
Limit of Detection (LoD)	0.11 ng/mL
Limit of Quantitation (LoQ)	0.26 ng/mL

### 15.2. Trueness

Trueness has been demonstrated through a recovery test for the Progesterone ELISA with the Progesterone standard from the National Metrology Institute of Japan (code 6003-a).

### 15.3. Precision

Precision of the Progesterone ELISA was determined by performing a complex precision study.

**Repeatability:** A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within run (Repeatability)	
			SD	CV%
1	75	0.93	0.08	9.0%
2	75	1.18	0.10	8.4%
3	75	2.60	0.22	8.5%
4	75	7.75	0.77	9.9%
5	75	22.05	1.60	7.2%
6	75	36.72	2.25	6.1%

**Reproducibility:** A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (ng/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	0.89	0.11	12.5%
2	150	1.14	0.15	13.0%
3	150	2.44	0.37	15.1%
4	150	7.41	1.06	14.3%
5	150	21.54	2.08	9.7%
6	150	36.43	2.33	6.4%

#### 15.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For progesterone concentration by Progesterone ELISA, the measurement procedure shows linearity for the interval from 0.17 to 41.02 ng/mL within the allowable deviation of linearity (ADL) of  $\pm 15\%$ .

#### 15.5. Method comparison

The Progesterone ELISA was compared against a commercially available quantitative ELISA, following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 109 samples, selected to represent a wide range of progesterone concentrations [0.17 – 35.84 ng/mL], was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
109	1.09 [0.99 to 1.21]	0.03 [-0.08 to 0.13]	0.98

#### 15.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested (unit)	Mean %Cross reactivity
Testosterone	1000 ng/mL	0.31%
17 $\alpha$ -OH Progesterone	1000 ng/mL	0.60%
17 $\beta$ -Estradiol	100 ng/mL	0.79%
Estrone	3000 ng/mL	0.00%
Estriol	1000 ng/mL	0.00%
Cortisol	100 $\mu$ g/mL	0.00%
Pregnenolone	100 ng/mL	1.71%
5 $\alpha$ -dihydrotestosterone	1000 ng/mL	0.32%
11-deoxycortisol	1000 ng/mL	0.05%
Corticosterone	1000 ng/mL	0.12%

The following substances do not interfere with a bias of  $> \pm 15\%$  in the Progesterone ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	10 g/dL
Triglycerides	500 mg/dL

#### 15.7. Serum-plasma study

The Progesterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 23 samples (all native) within the range of 0.26 to 40 ng/mL were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	0.93 [0.91 to 0.95]	0.01 [-0.20 to 0.21]	1.00
Lithium Heparin	0.95 [0.90 to 1.01]	0.16 [-0.41 to 0.74]	0.99
Sodium Heparin	0.95 [0.92 to 0.98]	0.12 [-0.20 to 0.45]	1.00
EDTA	0.98 [0.96 to 1.01]	-0.02 [-0.30 to 0.26]	1.00

## 16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays<sup>7</sup>. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

## 17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

## 18. BIBLIOGRAPHY

1. Wu SP, Li R, DeMayo FJ. Progesterone Receptor Regulation of Uterine Adaptation for Pregnancy. *Trends Endocrinol Metab.* 2018 Jul;29(7):481-491.
2. Kasum M, Simunić V, Vrčić H, Stanić P, Orešković S, Beketić-Orešković L. Follicular progesterone elevations with ovulation induction for IVF. *Gynecol Endocrinol.* 2014 Aug;30(8):537-41.
3. Di Renzo GC, Giardina I, Clerici G, Brillo E, Gerli S. Progesterone in normal and pathological

pregnancy. *Horm Mol Biol Clin Investig.* 2016 Jul 1;27(1):35-48.

4. Feng C, Chen ZY, Zhang J, Xu H, Zhang XM, Huang XF. Clinical utility of serum reproductive hormones for the early diagnosis of ectopic pregnancy in the first trimester. *J Obstet Gynaecol Res.* 2013 Feb;39(2):528-35.
5. Rausch ME, Sammel MD, Takacs P, Chung K, Shaunik A, Barnhart KT. Development of a multiple marker test for ectopic pregnancy. *Obstet Gynecol.* 2011 Mar;117(3):573-82.
6. Basic QC Practices On-line Course; <http://www.Westgard.com>.
7. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

## 19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

## 20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: [www.diametra.com](http://www.diametra.com).

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