DCM002-14

Ed. 02/2024

TESTOSTERONE ELISA

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18. BIBLIOGRAFIA

1. Jamerson. JL, de Kretser. D, Marshall. JC and De Groot. LJ. Endocrinology –adult and pediatric 6th edition. pp368-374.
2. Brambilla. DJ, Matsumoto. AM, Araujo. AB, McKinlay. JB. The Effect of Diurnal Variation on Clinical Measurement of Serum Testosterone and Other Sex Hormone Levels in Men. J Clin Endocrinol Metab. 2009 Mar; 94(3): 907–913.
3. Rajfer, J. Decreased Testosterone in the Aging Male. Rev Urol. 2003;5(suppl 1):S1–S2
4. Bhasin S, Brito JP, Cunningham GR, Hayes FJ, Hodis HN, Matsumoto AM, Snyder PJ, Swerdloff RS, Wu FC, Yialamas MA. Testosterone Therapy in Men With Hypogonadism: An Endocrine Society Clinical Practice Guideline. J Clin Endocrinol Metab. 2018 May 1;103(5):1715-1744.
5. Ho CK, Beckett GJ. Late-onset male hypogonadism: clinical and laboratory evaluation. J Clin Pathol. 2011 Jun;64(6):459-65.
6. Hugh Jones T. 'What should I do with a 60-year old man with a slightly low serum total testosterone concentration and normal levels of serum gonadotrophins'? Clin Endocrinol (Oxf). 2010 May;72(5):584-8.
7. Khera M, Adaikan G, Buvat J, Carrier S, El-Meliegy A, Hatzimouratidis K, McCullough A, Morgentaler A, Torres LO, Salonia A. Diagnosis and Treatment of Testosterone Deficiency: Recommendations From the Fourth International Consultation for Sexual Medicine (ICSM 2015). J Sex Med. 2016 Dec;13(12):1787-1804.
8. Dimopoulou C, Ceausu I, Depypere H, Lambrinoudaki I, Mueck A, Pérez-López FR, Rees M, van der Schouw YT, Senturk LM, Simonsini T, Stevenson JC, Stute P, Goulis DG. EMAS position statement: Testosterone replacement therapy in the aging male. Maturitas. 2016 Feb;84:94-9
9. Aversa A, Morgentaler A. The practical management of testosterone deficiency in men. Nat Rev Urol. 2015 Nov;12(11):641-50.
10. Basic QC Practices On-line Course; <http://www.Westgard.com>.
11. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27–33

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TESTOSTERONE ELISA

for routine analysis

Direct immunoenzymatic determination of total Testosterone in human serum and plasma

IVD

LOT

See external label

2°C  8°C $\Sigma = 96$ tests

REF DKO002

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use
For Laboratory Professional Use

Testosterone ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of total Testosterone in human serum or plasma from an adult population.

2. CLINICAL SIGNIFICANCE

Testosterone (17β -OH-4-androstene-3-one) has a molecular weight of 288 Daltons and is considered the principal androgen found in circulation of mature male mammals¹. In males it is synthesised and secreted by the Leydig cells located in the interstitium of the testis; in females, testosterone is produced in various locations such as the ovaries, adrenal gland and peripheral tissues, with additional production due to inter-conversion from other steroid hormones. Testosterone is found in circulation predominantly linked to carrier proteins, the most common of which being sex-hormone binding protein (SHBG). Only 1 – 2% of testosterone is biologically active, i.e. not bound to any plasmatic protein, and is referred to as free testosterone.

Testosterone plays a key part in the development of reproductive tissues and secondary sex characteristics in men. There is an observed and well documented circadian variation of testosterone levels in men with the circulating concentration being higher in the morning and declining throughout the day². Testosterone levels also decline in ageing males (andropause) and is often associated with loss of muscle and bone mass, leading to osteoporosis, loss of libido, erectile dysfunction, depression and impaired cognitive function³.

Disorders related to high levels or low levels of testosterone affect mainly biological processes related to gonads, sexual functionality and secondary sexual characteristics both in male and female.

Clinical conditions can be associated with either testosterone excess (hyperandrogenism), deficit (hypoandrogenism) or to alteration in the normal levels that occur in pre-pubertal and pubertal age (precocious puberty or delayed puberty). In all these disorders, quantification of

testosterone is used as an aid in diagnosis when considered in association with other clinical evidence and laboratory data.

Treatment of hypogonadism includes testosterone replacement therapy, treatment with aromatase inhibitors (which inhibit the conversion of testosterone to oestrogen) and treatment with human chorionic gonadotropin (which acts as an LH analogue stimulating endogenous testosterone secretion from testicular Leydig cells)⁴⁻⁹. These treatments should be monitored, via measurement of testosterone levels at pre-determined time intervals.

3. PRINCIPLE OF THE METHOD

The Testosterone ELISA is a competitive enzyme immunometric assay (ELISA) where testosterone (antigen) in the sample competes with the antigenic testosterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti- testosterone coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H_2SO_4) is added. The colour intensity is inversely proportional to the testosterone concentration in the sample.

Testosterone concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (5 vials, 1 mL each)

Human serum, ProClin >0.0015% and NaN_3 >0.1%

CAL0
CAL1
CAL2
CAL3
CAL4

REF	DCE002/0206-0
REF	DCE002/0207-0
REF	DCE002/0208-0
REF	DCE002/0209-0
REF	DCE002/0210-0

2. Control (2 vials, 1 mL each)

Human serum, ProClin >0.0015% and NaN_3 >0.1%

Control A

REF DCE045/0203A-0

Control B

REF DCE045/0203B-0

Controls Concentration is indicated on the Certificate of Analysis

3. Conjugate (1 vial, 12 mL)

Testosterone conjugated with horseradish peroxidase (HRP), contains Proclin > 0.0015%, BSA, sheep serum

REF DCE002/0202-0

4. Coated Microplate (1 breakable microplate)

Anti Testosterone antibody adsorbed on microplate

REF DCE002/0203-0

5. TMB Substrate (1 vial, 15 mL)

H_2O_2 -TMB 0.26 g/L (avoid any skin contact)

ProClin <0.0015%

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

REF DCE005-0

7. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4, contains Proclin > 0.0015%

REF DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, controls, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.

- Classification according to Regulation (EC) No. 1272/2008 [CLP]
Skin sensitivity, Category 1

Contains: ProClin 300

Warning

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing

dust/fume/gas/mist/vapours/spray.

P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- Some reagents (calibrators and controls) contain small amounts of Sodium Azide (NaN_3) <0.1% which may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through with large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of

samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- Once opened, the calibrators are stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	24 hours
Freeze/thaw cycles	1 cycle

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

The calibrators are ready for use and have the following concentration of testosterone:

	C ₀	C ₁	C ₂	C ₃	C ₄
ng/mL	0	0.6	1.4	4.0	16.0

The controls are ready to use; the concentrations of the controls are printed on the label.

9.2. Preparation of the Conjugate

Ready to use. Mix gently for 5 minutes on a roller mixer.

9.3. Preparation of the Wash Solution

Dilute the content of the vial “10X Conc. Wash Solution” with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of total testosterone can be performed in human serum or plasma samples.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation. Testing of heat-inactivated sera is not recommended.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Sample/ Control		25 µL	
Calibrator C ₀ -C ₄	25 µL		
Conjugate	100 µL	100 µL	

Incubate at 37°C for 1 hour.

Remove the content from each well. Wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: in case you use an automatic washer, it is advised to do 6 washing steps.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at room temperature (22 – 28°C) for 15 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories¹⁰.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required**. A smoothed spline fit including Calibrator 0 can be used. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

Conversion of units

To convert results to SI units:

$$\text{pmol/L} = \text{pg/mL} \times 3.47$$

To convert results to mass units:

$$\text{pg/mL} = \text{pmol/L} \times 0.288$$

12. MEASURING RANGE

The assay measuring range (AMR) is 0.30 – 16.0 ng/mL (1.04 – 55.47 nmol/L).

Any value that reads below 0.30 ng/mL (1.04 nmol/L) should be reported as "< 0.30 ng/mL (1.04 nmol/L)". Any value that reads above 16.0 ng/mL (55.47 nmol/L) should be reported as "> 16.0 ng/mL (55.47 nmol/L)".

13. METROLOGY AND TRACEABILITY

The calibrators of this kit are traceable to the National Metrology Institute of Australia (NMIA) standard for testosterone (code: M914c).

14. EXPECTED VALUES

The following ranges were determined using the Testosterone ELISA and are provided for information only. The 95% reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 "Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory".

Adults	n	Median (ng/mL)	Expected values (ng/mL)
Males			
21 – 49 years	120	4.75	2.45 – 8.79
> 50 years	120	4.14	1.50 – 6.92
Females			
Premenopausal	122	<LoQ	< LoQ – 0.56
Postmenopausal	120	0.40	< LoQ – 0.66

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.05 ng/mL
Limit of Detection (LoD)	0.17 ng/mL
Limit of Quantitation (LoQ)	0.30 ng/mL

15.2. Trueness

Trueness has been demonstrated through a recovery test for the Testosterone ELISA with the National Metrology Institute of Australia (NMI) standard for testosterone (code: M914c).

15.3. Precision

Precision of the Testosterone ELISA was determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within run (Repeatability)	
			SD	CV%
1	75	0.89	0.06	7.1%
2	75	1.42	0.10	7.4%
3	75	2.87	0.15	5.2%
4	75	4.97	0.20	3.9%
5	75	7.77	0.24	3.2%
6	75	14.71	0.53	3.6%

Reproducibility: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (ng/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	0.88	0.09	10.5%
2	150	1.40	0.19	13.7%
3	150	2.86	0.35	12.3%
4	150	5.01	0.47	9.3%
5	150	7.68	0.59	7.7%
6	150	14.19	1.01	7.1%

15.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For testosterone concentration by Testosterone ELISA, the measurement procedure shows linearity for the interval from 0.30 to 17.30 ng/mL (1.04 to 59.98 nmol/L) within the allowable deviation of linearity (ADL) of $\pm 15\%$.

15.5. Method comparison

The Testosterone ELISA was compared against a commercially available quantitative Testosterone ELISA, following CLSI EP-9C, "Method Comparison and Bias Estimation Using Patient Samples". A total of 52 samples, selected to represent a wide range of testosterone concentrations was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
52	1.06 [1.01 to 1.13]	0.33 [0.27 to 0.48]	0.99

15.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested	Mean %Cross reactivity
5 α -dihydrotestosterone	500 ng/mL	2.0%
Androstenedione	100 ng/mL	4.4%
Androsterone	100 ng/mL	0.1%
DHEA-S	10000 ng/mL	0.0%
Cortisol	1000 ng/mL	0.0%
Cortisone	1000 ng/mL	0.0%
17B-Estradiol	5 μ g/mL	0.0%
Prednisone	1000 ng/mL	0.0%
Estrone	1000 ng/mL	-0.1%
Ethisterone	100 ng/mL	4.3%
Testosterone propionate	1000 ng/mL	0.8%
11-keto-testosterone	100 ng/mL	10.0%
11-B-hydroxy-testosterone	10 ng/mL	69.7%
Prednisolone	200 ng/mL	0.0%
Progesterone	1000 ng/mL	0.0%

Danazol	1000 ng/mL	1.2%
Dexamethasone	2000 ng/mL	0.0%
Abiraterone	400 ng/mL	-0.2%

The following substances do not interfere with a bias of $> \pm 15\%$ in the Testosterone ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	7 g/dL
Triglyceride	500 mg/dL

15.7. Serum-plasma study

The Testosterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP35-Ed1 guidelines. A total of 20 samples (16 native, 4 spiked) to cover the assay range of were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	0.98 [0.93 to 1.02]	0.10 [-0.19 to 0.39]	1.00
Lithium Heparin	1.01 [0.97 to 1.05]	0.03 [-0.21 to 0.26]	1.00
Sodium Heparin	1.00 [0.97 to 1.04]	0.05 [-0.17 to 0.28]	1.00
EDTA	1.01 [0.96 to 1.07]	0.19 [-0.15 to 0.53]	0.99

16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays¹¹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

18. BIBLIOGRAPHY

1. Jamerson. JL, de Kretser. D, Marshall. JC and De Groot. LJ. Endocrinology –adult and pediatric 6th edition. pp368-374.
2. Brambilla. DJ, Matsumoto. AM, Araujo. AB, McKinlay. JB. The Effect of Diurnal Variation on Clinical Measurement of Serum Testosterone and Other Sex Hormone Levels in Men. J Clin Endocrinol Metab. 2009 Mar; 94(3): 907–913.
3. Rajfer, J. Decreased Testosterone in the Aging Male. Rev Urol. 2003;5(suppl 1):S1–S2
4. Bhasin S, Brito JP, Cunningham GR, Hayes FJ, Hodis HN, Matsumoto AM, Snyder PJ, Swerdloff RS, Wu FC, Yialamas MA. Testosterone Therapy in Men With Hypogonadism: An Endocrine Society Clinical Practice Guideline. J Clin Endocrinol Metab. 2018 May 1;103(5):1715-1744.
5. Ho CK, Beckett GJ. Late-onset male hypogonadism: clinical and laboratory evaluation. J Clin Pathol. 2011 Jun;64(6):459-65.
6. Hugh Jones T. 'What should I do with a 60-year old man with a slightly low serum total testosterone concentration and normal levels of serum gonadotrophins'? Clin Endocrinol (Oxf). 2010 May;72(5):584-8.
7. Khera M, Adaikan G, Buvat J, Carrier S, El-Meliegy A, Hatzimouratidis K, McCullough A, Morgentaler A, Torres LO, Salonia A. Diagnosis and Treatment of Testosterone Deficiency: Recommendations From the Fourth International Consultation for Sexual Medicine (ICSM 2015). J Sex Med. 2016 Dec;13(12):1787-1804.
8. Dimopoulou C, Ceausu I, Depypere H, Lambrinoudaki I, Mueck A, Pérez-López FR, Rees M, van der Schouw YT, Senturk LM, Simonsini T, Stevenson JC, Stute P, Goulis DG. EMAS position statement: Testosterone replacement therapy in the aging male. Maturitas. 2016 Feb;84:94-9
9. Aversa A, Morgentaler A. The practical management of testosterone deficiency in men. Nat Rev Urol. 2015 Nov;12(11):641-50.
10. Basic QC Practices On-line Course; <http://www.Westgard.com>.
11. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27–33

19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

Ed. 02/2024

DCM002-14

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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